



LIPOSUCTION

CANNULAS & ACCESSORIES



LOTUSMED
INSTRUMENTS

www.lotusmedinstruments.com
Email: lotusmedinstruments@gmail.com
Whatsapp & Call: +92 317 7290854



Table of Contents

- A. Liposuction Cannulas
 - a) Infiltration Cannulas
 - b) General Suction Cannulas
 - c) Harvesting Cannulas
 - e) Injecting Cannulas
 - f) Speciality Cannulas

- B. Accessories

- C. Our Specialized Sets

Cannula Hubs/ Connections



Luer lock



Luer lock
(Stainless Steel)



Super Luer Lock



Treaded



60cc Tommey hub

Machine Connections



ML-1000.01

Cannulas and Harvesters as a complete instrument ready to use with an aspiration machine for liposuction and large volume fat harvesting procedures.

Fixed Handle

Luer Lock Hub/Connector:

The universal Luer Lock hub enables cannula to attach directly with any Luer Lock Syringe.

Super Luer Lock Hub/Connector:

A modification to general Luer Lock hub which is more powerful, more stable, Leak proof, free from instability and contamination.

Threaded Hub/Connector:

Threaded connector hub style is used mostly for machine assisted Liposuction cannulas, Medicta's unique Light weight threaded cannula can also be used with ordinary syringe. All you need is Medicta's connector from threaded to Luer style.


60cc Tommey Hub/Connector:

An ideal hub style for surgeons why love to handle large fat volumes. This hub style connect cannula to a 60cc Tommey to syringe.

Fixed Handle Cannulas:

The unique design of handle offer powerful grip harvesting large volumes. You will feel amazingly comfortable ever after a long operational time.


The thumb grip provide a unique opportunity to control the thrust.



INFILTRATION CANNULAS

Infiltration Cannulas are designed for safe and effective delivery of tumescent solution. The spiral micro holes on cannula ensure an even distribution of local anesthetic which results in a longer anesthetic effect, this is less painful and more pleasant liposuction experience for patient.

INFILTRATION CANNULAS



Tumescent Infiltrator

This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

60cc Tommey hub

LT-101.01	Ø 2.5 mm X 15 cm
LT-101.02	Ø 2.5 mm X 20 cm
LT-101.03	Ø 3.0 mm X 20 cm
LT-101.04	Ø 3.0 mm X 25 cm

Treaded Hub

LT-102.01	Ø 2.0 mm X 15 cm
LT-102.02	Ø 3.0 mm X 15 cm
LT-102.03	Ø 3.0 mm X 20 cm
LT-102.04	Ø 3.0 mm X 30 cm
LT-102.05	Ø 4.0 mm X 20 cm
LT-102.06	Ø 4.0 mm X 30 cm

Super Luer Lock

LT-103.01	Ø 1.0 mm X 5 cm
LT-103.02	Ø 1.0 mm X 7 cm
LT-103.03	Ø 1.5 mm X 7 cm
LT-103.04	Ø 1.5 mm X 9 cm

Luer lock

LT-104.01	Ø 1.0 mm X 4 cm
LT-104.02	Ø 1.0 mm X 5 cm
LT-104.03	Ø 1.5 mm X 7 cm
LT-104.04	Ø 1.5 mm X 9 cm
LT-104.05	Ø 1.5 mm X 10 cm
LT-104.06	Ø 1.5 mm X 12 cm
LT-104.07	Ø 2.0 mm X 10 cm
LT-104.08	Ø 2.0 mm X 12 cm
LT-104.09	Ø 2.0 mm X 15 cm
LT-104.10	Ø 2.0 mm X 20 cm
LT-104.11	Ø 2.5 mm X 15 cm
LT-104.12	Ø 2.5 mm X 20 cm
LT-104.13	Ø 3.0 mm X 15 cm
LT-104.14	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)

LT-105.01	Ø 1.0 mm X 4 cm
LT-105.02	Ø 1.0 mm X 5 cm
LT-105.03	Ø 1.0 mm X 7 cm
LT-105.04	Ø 1.5 mm X 7 cm
LT-105.05	Ø 1.5 mm X 9 cm
LT-105.06	Ø 1.5 mm X 10 cm
LT-105.07	Ø 1.5 mm X 12 cm
LT-105.08	Ø 2.0 mm X 10 cm
LT-105.09	Ø 2.0 mm X 12 cm
LT-105.10	Ø 2.0 mm X 15 cm
LT-105.11	Ø 2.0 mm X 20 cm
LT-105.12	Ø 2.5 mm X 15 cm
LT-105.13	Ø 2.5 mm X 20 cm
LT-105.14	Ø 3.0 mm X 15 cm
LT-105.15	Ø 3.0 mm X 20 cm

Power Cannula

LT-106.01	Ø 2.0 mm X 15 cm
LT-106.02	Ø 2.0 mm X 20 cm
LT-106.03	Ø 2.0 mm X 25 cm
LT-106.04	Ø 2.0 mm X 30 cm
LT-106.05	Ø 2.0 mm X 35 cm
LT-106.06	Ø 3.0 mm X 15 cm
LT-106.07	Ø 3.0 mm X 20 cm
LT-106.08	Ø 3.0 mm X 25 cm
LT-106.09	Ø 3.0 mm X 30 cm
LT-106.10	Ø 3.0 mm X 35 cm
LT-106.11	Ø 4.0 mm X 15 cm
LT-106.12	Ø 4.0 mm X 20 cm
LT-106.13	Ø 4.0 mm X 25 cm
LT-106.14	Ø 4.0 mm X 30 cm
LT-106.15	Ø 4.0 mm X 35 cm
LT-106.16	Ø 5.0 mm X 15 cm
LT-106.17	Ø 5.0 mm X 20 cm
LT-106.18	Ø 5.0 mm X 25 cm
LT-106.19	Ø 5.0 mm X 30 cm
LT-106.20	Ø 5.0 mm X 35 cm
LT-106.21	Ø 6.0 mm X 15 cm
LT-106.22	Ø 6.0 mm X 20 cm
LT-106.23	Ø 6.0 mm X 25 cm
LT-106.24	Ø 6.0 mm X 30 cm
LT-106.25	Ø 6.0 mm X 35 cm

INFILTRATION CANNULAS



Tumescent Bayonet Infiltrator

The bayonet shaped handles allow an easier manipulation

60cc Tommey hub

LT-101.11	Ø 2.0 mm X 15 cm
LT-101.12	Ø 2.5 mm X 20 cm
LT-101.13	Ø 3.0 mm X 20 cm
LT-101.14	Ø 3.0 mm X 25 cm

Treaded Hub

LT-102.11	Ø 2.0mm X 15 cm
LT-102.12	Ø 3.0mm X 15 cm
LT-102.13	Ø 3.0mm X 20 cm
LT-102.14	Ø 3.0mm X 30 cm
LT-102.15	Ø 4.0mm X 20 cm
LT-102.16	Ø 4.0mm X 30 cm

Super Luer Lock

LT-103.11	Ø 1.0 mm X 5 cm
LT-103.12	Ø 1.0 mm X 7 cm
LT-103.13	Ø 1.5 mm X 7 cm
LT-103.14	Ø 1.5 mm X 9 cm

Luer lock

LT-104.11	Ø 1.0 mm X 4 cm
LT-104.12	Ø 1.0 mm X 5 cm
LT-104.13	Ø 1.0 mm X 7 cm
LT-104.14	Ø 1.5 mm X 7 cm
LT-104.15	Ø 1.5 mm X 9 cm
LT-104.16	Ø 1.5 mm X 10 cm
LT-104.17	Ø 1.5 mm X 12 cm
LT-104.18	Ø 2.0 mm X 10 cm
LT-104.19	Ø 2.0 mm X 12 cm
LT-104.20	Ø 2.0 mm X 15 cm
LT-104.21	Ø 2.0 mm X 20 cm
LT-104.22	Ø 2.5 mm X 15 cm
LT-104.23	Ø 2.5 mm X 20 cm
LT-104.24	Ø 3.0 mm X 15 cm
LT-104.25	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)

LT-105.21	Ø 1.0 mm X 4 cm
LT-105.22	Ø 1.0 mm X 5 cm
LT-105.23	Ø 1.0 mm X 7 cm
LT-105.24	Ø 1.5 mm X 7 cm
LT-105.25	Ø 1.5 mm X 9 cm
LT-105.26	Ø 1.5 mm X 10 cm
LT-105.27	Ø 1.5 mm X 12 cm
LT-105.28	Ø 2.0 mm X 10 cm
LT-105.29	Ø 2.0 mm X 12 cm
LT-105.30	Ø 2.0 mm X 15 cm
LT-105.31	Ø 2.0 mm X 20 cm
LT-105.32	Ø 2.5 mm X 15 cm
LT-105.33	Ø 2.5 mm X 20 cm
LT-105.34	Ø 3.0 mm X 15 cm
LT-105.35	Ø 3.0 mm X 20 cm

Power Cannula

LT-106.31	Ø 2 mm X 15 cm
LT-106.32	Ø 2 mm X 20 cm
LT-106.33	Ø 2 mm X 25 cm
LT-106.34	Ø 2 mm X 30 cm
LT-106.35	Ø 2 mm X 35 cm
LT-106.36	Ø 3 mm X 15 cm
LT-106.37	Ø 3 mm X 20 cm
LT-106.38	Ø 3 mm X 25 cm
LT-106.39	Ø 3 mm X 30 cm
LT-106.40	Ø 3 mm X 35 cm
LT-106.41	Ø 4 mm X 15 cm
LT-106.42	Ø 4 mm X 20 cm
LT-106.43	Ø 4 mm X 25 cm
LT-106.44	Ø 4 mm X 30 cm
LT-106.45	Ø 4 mm X 35 cm
LT-106.46	Ø 5 mm X 15 cm
LT-106.47	Ø 5 mm X 20 cm
LT-106.48	Ø 5 mm X 25 cm
LT-106.49	Ø 5 mm X 30 cm
LT-106.50	Ø 5 mm X 35 cm
LT-106.51	Ø 6 mm X 15 cm
LT-106.52	Ø 6 mm X 20 cm
LT-106.53	Ø 6 mm X 25 cm
LT-106.54	Ø 6 mm X 30 cm
LT-106.55	Ø 6 mm X 35 cm

INFILTRATION CANNULAS



Tumescent Infiltrator Angled

This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.

60cc Tommey hub

LT-101.21	Ø 2.0 mm X 15 cm
LT-101.22	Ø 2.5 mm X 20 cm
LT-101.23	Ø 3.0 mm X 20 cm
LT-101.24	Ø 3.0 mm X 25 cm

Treaded Hub

LT-102.21	Ø 2.0mm X 15 cm
LT-102.22	Ø 3.0mm X 15 cm
LT-102.23	Ø 3.0mm X 20 cm
LT-102.24	Ø 3.0mm X 30 cm
LT-102.25	Ø 4.0mm X 20 cm
LT-102.26	Ø 4.0mm X 30 cm

Super Luer Lock

LT-103.21	Ø 1.0 mm X 5 cm
LT-103.22	Ø 1.0 mm X 7 cm
LT-103.23	Ø 1.5 mm X 7 cm
LT-103.24	Ø 1.5 mm X 9 cm

Luer lock

LT-104.31	Ø 1.0 mm X 4 cm
LT-104.32	Ø 1.0 mm X 5 cm
LT-104.33	Ø 1.0 mm X 7 cm
LT-104.34	Ø 1.5 mm X 7 cm
LT-104.35	Ø 1.5 mm X 9 cm
LT-104.36	Ø 1.5 mm X 10 cm
LT-104.37	Ø 1.5 mm X 12 cm
LT-104.38	Ø 2.0 mm X 10 cm
LT-104.39	Ø 2.0 mm X 12 cm
LT-104.40	Ø 2.0 mm X 15 cm
LT-104.41	Ø 2.0 mm X 20 cm
LT-104.42	Ø 2.5 mm X 15 cm
LT-104.43	Ø 2.5 mm X 20 cm
LT-104.44	Ø 3.0 mm X 15 cm
LT-104.45	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)

LT-105.41	Ø 1.0 mm X 4 cm
LT-105.42	Ø 1.0 mm X 5 cm
LT-105.43	Ø 1.0 mm X 7 cm
LT-105.44	Ø 1.5 mm X 7 cm
LT-105.45	Ø 1.5 mm X 9 cm
LT-105.46	Ø 1.5 mm X 10 cm
LT-105.47	Ø 1.5 mm X 12 cm
LT-105.48	Ø 2.0 mm X 10 cm
LT-105.49	Ø 2.0 mm X 12 cm
LT-105.50	Ø 2.0 mm X 15 cm
LT-105.51	Ø 2.0 mm X 20 cm
LT-105.52	Ø 2.5 mm X 15 cm
LT-105.53	Ø 2.5 mm X 20 cm
LT-105.54	Ø 3.0 mm X 15 cm
LT-105.55	Ø 3.0 mm X 20 cm

Power Cannula

LT-106.61	Ø 2 mm X 15 cm
LT-106.62	Ø 2 mm X 20 cm
LT-106.63	Ø 2 mm X 25 cm
LT-106.64	Ø 2 mm X 30 cm
LT-106.65	Ø 2 mm X 35 cm
LT-106.66	Ø 3 mm X 15 cm
LT-106.67	Ø 3 mm X 20 cm
LT-106.68	Ø 3 mm X 25 cm
LT-106.69	Ø 3 mm X 30 cm
LT-106.70	Ø 3 mm X 35 cm
LT-106.71	Ø 4 mm X 15 cm
LT-106.72	Ø 4 mm X 20 cm
LT-106.73	Ø 4 mm X 25 cm
LT-106.74	Ø 4 mm X 30 cm
LT-106.75	Ø 4 mm X 35 cm
LT-106.76	Ø 5 mm X 15 cm
LT-106.77	Ø 5 mm X 20 cm
LT-106.78	Ø 5 mm X 25 cm
LT-106.79	Ø 5 mm X 30 cm
LT-106.80	Ø 5 mm X 35 cm
LT-106.81	Ø 6 mm X 15 cm
LT-106.82	Ø 6 mm X 20 cm
LT-106.83	Ø 6 mm X 25 cm
LT-106.84	Ø 6 mm X 30 cm
LT-106.85	Ø 6 mm X 35 cm

INFILTRATION CANNULAS



Facelift Infiltration Cannula

This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.



60cc Tommey hub	
LT-101.31	Ø 2.0 mm X 15 cm
LT-101.32	Ø 2.5 mm X 20 cm
LT-101.33	Ø 3.0 mm X 20 cm
LT-101.34	Ø 3.0 mm X 25 cm


Treaded Hub	
LT-102.31	Ø 2.0mm X 15 cm
LT-102.32	Ø 3.0mm X 15 cm
LT-102.33	Ø 3.0mm X 20 cm
LT-102.34	Ø 3.0mm X 30 cm
LT-102.35	Ø 4.0mm X 20 cm
LT-102.36	Ø 4.0mm X 30 cm

Super Luer Lock	
LT-103.31	Ø 1.0 mm X 5 cm
LT-103.32	Ø 1.0 mm X 7 cm
LT-103.33	Ø 1.5 mm X 7 cm
LT-103.34	Ø 1.5 mm X 9 cm

Luer lock	
LT-104.51	Ø 1.0 mm X 4 cm
LT-104.52	Ø 1.0 mm X 5 cm
LT-104.53	Ø 1.0 mm X 7 cm
LT-104.54	Ø 1.5 mm X 7 cm
LT-104.55	Ø 1.5 mm X 9 cm
LT-104.56	Ø 1.5 mm X 10 cm
LT-104.57	Ø 1.5 mm X 12 cm
LT-104.58	Ø 2.0 mm X 10 cm
LT-104.59	Ø 2.0 mm X 12 cm
LT-104.50	Ø 2.0 mm X 15 cm
LT-104.51	Ø 2.0 mm X 20 cm
LT-104.52	Ø 2.5 mm X 15 cm
LT-104.53	Ø 2.5 mm X 20 cm
LT-104.54	Ø 3.0 mm X 15 cm
LT-104.55	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)	
LT-105.61	Ø 1.0 mm X 4 cm
LT-105.62	Ø 1.0 mm X 5 cm
LT-105.63	Ø 1.0 mm X 7 cm
LT-105.64	Ø 1.5 mm X 7 cm
LT-105.65	Ø 1.5 mm X 9 cm
LT-105.66	Ø 1.5 mm X 10 cm
LT-105.67	Ø 1.5 mm X 12 cm
LT-105.68	Ø 2.0 mm X 10 cm
LT-105.69	Ø 2.0 mm X 12 cm
LT-105.70	Ø 2.0 mm X 15 cm
LT-105.71	Ø 2.0 mm X 20 cm
LT-105.72	Ø 2.5 mm X 15 cm
LT-105.73	Ø 2.5 mm X 20 cm
LT-105.74	Ø 3.0 mm X 15 cm
LT-105.75	Ø 3.0 mm X 20 cm

Power Cannula	
LT-106.91	Ø 2 mm X 15 cm
LT-106.92	Ø 2 mm X 20 cm
LT-106.93	Ø 2 mm X 25 cm
LT-106.94	Ø 2 mm X 30 cm
LT-106.95	Ø 2 mm X 35 cm
LT-106.96	Ø 3 mm X 15 cm
LT-106.97	Ø 3 mm X 20 cm
LT-106.98	Ø 3 mm X 25 cm
LT-106.99	Ø 3 mm X 30 cm
LT-106.100	Ø 3 mm X 35 cm
LT-106.101	Ø 4 mm X 15 cm
LT-106.102	Ø 4 mm X 20 cm
LT-106.103	Ø 4 mm X 25 cm
LT-106.104	Ø 4 mm X 30 cm
LT-106.105	Ø 4 mm X 35 cm
LT-106.106	Ø 5 mm X 15 cm
LT-106.107	Ø 5 mm X 20 cm
LT-106.108	Ø 5 mm X 25 cm
LT-106.109	Ø 5 mm X 30 cm
LT-106.110	Ø 5 mm X 35 cm
LT-106.111	Ø 6 mm X 15 cm
LT-106.112	Ø 6 mm X 20 cm
LT-106.113	Ø 6 mm X 25 cm
LT-106.114	Ø 6 mm X 30 cm
LT-106.115	Ø 6 mm X 35 cm



GENERAL SUCTION CANNULA

General suction cannula are ideal tool for large volume Harvesting with ease perfection. These cannulas are designed to facilitate you to get perfectly designed goals in minimum time.

GENERAL SUCTION CANNULAS



Stevens Style Cannula

60cc Tommey hub	
LT-201.01	Ø 2.0 mm X 15 cm
LT-201.02	Ø 2.5 mm X 20 cm
LT-201.03	Ø 3.0 mm X 20 cm
LT-201.04	Ø 3.0 mm X 25 cm

Treaded Hub	
LT-202.01	Ø 3.0 mm X 20 cm
LT-202.02	Ø 3.0 mm X 30 cm
LT-202.03	Ø 4.0 mm X 30 cm
LT-202.04	Ø 5.0 mm X 30 cm
LT-202.05	Ø 6.0 mm X 30 cm


Super Luer Lock	
LT-203.01	Ø 1.0 mm X 5 cm
LT-203.02	Ø 1.0 mm X 7 cm
LT-203.03	Ø 1.5 mm X 7 cm
LT-203.04	Ø 1.5 mm X 9 cm
LT-203.05	Ø 1.5 mm X 7 cm

Luer lock	
LT-204.01	Ø 1.0 mm X 4 cm
LT-204.02	Ø 1.0 mm X 5 cm
LT-204.03	Ø 1.0 mm X 7 cm
LT-204.04	Ø 1.5 mm X 7 cm
LT-204.05	Ø 1.5 mm X 9 cm
LT-204.06	Ø 1.5 mm X 10 cm
LT-204.07	Ø 1.5 mm X 12 cm
LT-204.08	Ø 2.0 mm X 10 cm
LT-204.09	Ø 2.0 mm X 12 cm
LT-204.10	Ø 2.0 mm X 15 cm
LT-204.11	Ø 2.0 mm X 20 cm
LT-204.12	Ø 2.5 mm X 15 cm
LT-204.13	Ø 2.5 mm X 20 cm
LT-204.14	Ø 3.0 mm X 15 cm
LT-204.15	Ø 3.0 mm X 20 cm

Luer lock (Stainless Steel)	
LT-205.01	Ø 1.0 mm X 4 cm
LT-205.02	Ø 1.0 mm X 5 cm
LT-205.03	Ø 1.0 mm X 7 cm
LT-205.04	Ø 1.5 mm X 7 cm
LT-205.05	Ø 1.5 mm X 9 cm
LT-205.06	Ø 1.5 mm X 10 cm
LT-205.07	Ø 1.5 mm X 12 cm
LT-205.08	Ø 2.0 mm X 10 cm
LT-205.09	Ø 2.0 mm X 12 cm
LT-205.10	Ø 2.0 mm X 15 cm
LT-205.11	Ø 2.0 mm X 20 cm
LT-205.12	Ø 2.5 mm X 15 cm
LT-205.13	Ø 2.5 mm X 20 cm
LT-205.14	Ø 3.0 mm X 15 cm
LT-205.15	Ø 3.0 mm X 20 cm

Power Cannula	
LT-206.01	Ø 2 mm X 15 cm
LT-206.02	Ø 2 mm X 20 cm
LT-206.03	Ø 2 mm X 25 cm
LT-206.04	Ø 2 mm X 30 cm
LT-206.05	Ø 2 mm X 35 cm
LT-206.06	Ø 3 mm X 15 cm
LT-206.07	Ø 3 mm X 20 cm
LT-206.08	Ø 3 mm X 25 cm
LT-206.09	Ø 3 mm X 30 cm
LT-206.10	Ø 3 mm X 35 cm
LT-206.11	Ø 4 mm X 15 cm
LT-206.12	Ø 4 mm X 20 cm
LT-206.13	Ø 4 mm X 25 cm
LT-206.14	Ø 4 mm X 30 cm
LT-206.15	Ø 4 mm X 35 cm
LT-206.16	Ø 5 mm X 15 cm
LT-206.17	Ø 5 mm X 20 cm
LT-206.18	Ø 5 mm X 25 cm
LT-206.19	Ø 5 mm X 30 cm
LT-206.20	Ø 5 mm X 35 cm
LT-206.21	Ø 6 mm X 15 cm
LT-206.22	Ø 6 mm X 20 cm
LT-206.23	Ø 6 mm X 25 cm
LT-206.24	Ø 6 mm X 30 cm
LT-206.25	Ø 6 mm X 35 cm

GENERAL SUCTION CANNULAS



Cobra Bibevel Style

More aggressive than the Standard, Pyramid or Mercedes, this cannula has a tapered tip with two ports and works effectively within dense fibrous tissue.

60cc Tommey hub

LT-201.11	Ø 2.5 mm X 15 cm5
LT-201.12	Ø 2.5 mm X 20 cm
LT-201.13	Ø 3.0 mm X 15 cm
LT-201.14	Ø 3.0 mm X 20 cm
LT-201.15	Ø 3.0 mm X 25 cm
LT-201.16	Ø 4.0 mm X 25 cm
LT-201.17	Ø 4.0 mm X 30 cm
LT-201.18	Ø 4.0 mm X 35 cm
LT-201.19	Ø 4.5 mm X 30 cm
LT-201.20	Ø 4.5 mm X 35 cm

Treaded Hub

LT-202.11	Ø 3.0 mm X 15 cm
LT-202.12	Ø 3.0 mm X 25 cm
LT-202.13	Ø 4.0 mm X 15 cm
LT-202.14	Ø 4.0 mm X 25 cm

Super Luer Lock

LT-203.11	Ø 2.0 mm X 15 cm
LT-203.12	Ø 2.5 mm X 15 cm
LT-203.13	Ø 2.5 mm X 20 cm
LT-203.14	Ø 3.0 mm X 15 cm
LT-203.15	Ø 3.0 mm X 20 cm
LT-203.16	Ø 3.0 mm X 25 cm

Luer lock

LT-204.21	Ø 2.0 mm X 15 cm
LT-204.22	Ø 2.5 mm X 15 cm
LT-204.23	Ø 2.5 mm X 20 cm
LT-204.24	Ø 3.0 mm X 15 cm
LT-204.25	Ø 3.0 mm X 20 cm
LT-204.26	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

LT-205.21	Ø 2.0 mm X 15 cm
LT-205.22	Ø 2.5 mm X 15 cm
LT-205.23	Ø 2.5 mm X 20 cm
LT-205.24	Ø 3.0 mm X 15 cm
LT-205.25	Ø 3.0 mm X 20 cm
LT-205.26	Ø 3.0 mm X 25 cm

Power Cannula

LT-206.31	Ø 2mm X 15 cm
LT-206.32	Ø 2mm X 20 cm
LT-206.33	Ø 2mm X 25 cm
LT-206.34	Ø 2mm X 30 cm
LT-206.35	Ø 2mm X 35 cm
LT-206.36	Ø 3mm X 15 cm
LT-206.37	Ø 3mm X 20 cm
LT-206.38	Ø 3mm X 25 cm
LT-206.39	Ø 3mm X 30 cm
LT-206.40	Ø 3mm X 35 cm
LT-206.41	Ø 4mm X 15 cm
LT-206.42	Ø 4mm X 20 cm
LT-206.43	Ø 4mm X 25 cm
LT-206.44	Ø 4mm X 30 cm
LT-206.45	Ø 4mm X 35 cm
LT-206.46	Ø 5mm X 15 cm
LT-206.47	Ø 5mm X 20 cm
LT-206.48	Ø 5mm X 25 cm
LT-206.49	Ø 5mm X 30 cm
LT-206.50	Ø 5mm X 35 cm
LT-206.51	Ø 6mm X 15 cm
LT-206.52	Ø 6mm X 20 cm
LT-206.53	Ø 6mm X 25 cm
LT-206.54	Ø 6mm X 30 cm
LT-206.55	Ø 6mm X 35 cm

GENERAL SUCTION CANNULAS



Cobra Style Round Tip

Less aggressive than the Cobra Bibevel this cell friendly cannula has a rounded tip with two ports and works effectively within dense fibrous tissue during liposuction procedures.

60cc Tommey hub	
LT-201.21	Ø 2.5 mm X 15 cm
LT-201.22	Ø 2.5 mm X 20 cm
LT-201.23	Ø 3.0 mm X 15 cm
LT-201.24	Ø 3.0 mm X 20 cm
LT-201.25	Ø 3.0 mm X 25 cm
LT-201.26	Ø 4.0 mm X 25 cm
LT-201.27	Ø 4.0 mm X 30 cm
LT-201.28	Ø 4.0 mm X 35 cm
LT-201.29	Ø 4.5 mm X 30 cm
LT-201.30	Ø 4.5 mm X 35 cm

Treaded Hub	
LT-202.21	Ø 3.0 mm X 15 cm
LT-202.22	Ø 3.0 mm X 25 cm
LT-202.23	Ø 4.0 mm X 15 cm
LT-202.24	Ø 4.0 mm X 25 cm

Super Luer Lock	
LT-203.21	Ø 2.0 mm X 15 cm
LT-203.22	Ø 2.5 mm X 15 cm
LT-203.23	Ø 2.5 mm X 20 cm
LT-203.24	Ø 3.0 mm X 15 cm
LT-203.25	Ø 3.0 mm X 20 cm
LT-203.26	Ø 3.0 mm X 25 cm

Luer lock	
LT-204.31	Ø 2.0 mm X 15 cm
LT-204.32	Ø 2.5 mm X 15 cm
LT-204.33	Ø 2.5 mm X 20 cm
LT-204.34	Ø 3.0 mm X 15 cm
LT-204.35	Ø 3.0 mm X 20 cm
LT-204.36	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)	
LT-205.31	Ø 2.0 mm X 15 cm
LT-205.32	Ø 2.5 mm X 15 cm
LT-205.33	Ø 2.5 mm X 20 cm
LT-205.34	Ø 3.0 mm X 15 cm
LT-205.35	Ø 3.0 mm X 20 cm
LT-205.36	Ø 3.0 mm X 25 cm

Power Cannula	
LT-206.61	Ø 2 mm X 15 cm
LT-206.62	Ø 2 mm X 20 cm
LT-206.63	Ø 2 mm X 25 cm
LT-206.64	Ø 2 mm X 30 cm
LT-206.65	Ø 2 mm X 35 cm
LT-206.66	Ø 3 mm X 15 cm
LT-206.67	Ø 3 mm X 20 cm
LT-206.68	Ø 3 mm X 25 cm
LT-206.69	Ø 3 mm X 30 cm
LT-206.70	Ø 3 mm X 35 cm
LT-206.71	Ø 4 mm X 15 cm
LT-206.72	Ø 4 mm X 20 cm
LT-206.73	Ø 4 mm X 25 cm
LT-206.74	Ø 4 mm X 30 cm
LT-206.75	Ø 4 mm X 35 cm
LT-206.76	Ø 5 mm X 15 cm
LT-206.77	Ø 5 mm X 20 cm
LT-206.78	Ø 5 mm X 25 cm
LT-206.79	Ø 5 mm X 30 cm
LT-206.80	Ø 5 mm X 35 cm
LT-206.81	Ø 6 mm X 15 cm
LT-206.82	Ø 6 mm X 20 cm
LT-206.83	Ø 6 mm X 25 cm
LT-206.84	Ø 6 mm X 30 cm
LT-206.85	Ø 6 mm X 35 cm

GENERAL SUCTION CANNULAS



Standard, Single Hole

This is the most common cell friendly cannula design for soft tissue removal. The single port is located 180 degrees from the thumb rest on the hub.

60cc Tommey hub

LT-201.31	Ø 2.5 mm X 15 cm
LT-201.32	Ø 2.5 mm X 20 cm
LT-201.33	Ø 3.0 mm X 15 cm
LT-201.34	Ø 3.0 mm X 20 cm
LT-201.35	Ø 3.0 mm X 25 cm
LT-201.36	Ø 4.0 mm X 25 cm
LT-201.37	Ø 4.0 mm X 30 cm
LT-201.38	Ø 4.0 mm X 35 cm
LT-201.39	Ø 4.5 mm X 30 cm
LT-201.40	Ø 4.5 mm X 35 cm

Treaded Hub

LT-202.31	Ø 3.0 mm X 15 cm
LT-202.32	Ø 3.0 mm X 25 cm
LT-202.33	Ø 4.0 mm X 15 cm
LT-202.34	Ø 4.0 mm X 25 cm

Super Luer Lock

LT-203.31	Ø 1.5 mm X 12 cm
LT-203.32	Ø 2.0 mm X 15 cm
LT-203.33	Ø 2.0 mm X 20 cm
LT-203.34	Ø 2.5 mm X 15 cm
LT-203.35	Ø 2.5 mm X 20 cm
LT-203.36	Ø 3.0 mm X 15 cm
LT-203.37	Ø 3.0 mm X 20 cm
LT-203.38	Ø 3.0 mm X 25 cm

Luer lock

LT-204.41	Ø 1.5 mm X 12 cm
LT-204.42	Ø 2.0 mm X 15 cm
LT-204.43	Ø 2.0 mm X 20 cm
LT-204.44	Ø 2.5 mm X 15 cm
LT-204.45	Ø 2.5 mm X 20 cm
LT-204.46	Ø 3.0 mm X 15 cm
LT-204.47	Ø 3.0 mm X 20 cm
LT-204.48	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

LT-205.41	Ø 1.5 mm X 12 cm
LT-205.42	Ø 2.0 mm X 15 cm
LT-205.43	Ø 2.0 mm X 20 cm
LT-205.44	Ø 2.5 mm X 15 cm
LT-205.45	Ø 2.5 mm X 20 cm
LT-205.46	Ø 3.0 mm X 15 cm
LT-205.47	Ø 3.0 mm X 20 cm
LT-205.48	Ø 3.0 mm X 25 cm

Power Cannula

LT-206.91	Ø 2 mm X 15 cm
LT-206.92	Ø 2 mm X 20 cm
LT-206.93	Ø 2 mm X 25 cm
LT-206.94	Ø 2 mm X 30 cm
LT-206.95	Ø 2 mm X 35 cm
LT-206.96	Ø 3 mm X 15 cm
LT-206.97	Ø 3 mm X 20 cm
LT-206.98	Ø 3 mm X 25 cm
LT-206.99	Ø 3 mm X 30 cm
LT-206.100	Ø 3 mm X 35 cm
LT-206.101	Ø 4 mm X 15 cm
LT-206.102	Ø 4 mm X 20 cm
LT-206.103	Ø 4 mm X 25 cm
LT-206.104	Ø 4 mm X 30 cm
LT-206.105	Ø 4 mm X 35 cm
LT-206.106	Ø 5 mm X 15 cm
LT-206.107	Ø 5 mm X 20 cm
LT-206.108	Ø 5 mm X 25 cm
LT-206.109	Ø 5 mm X 30 cm
LT-206.110	Ø 5 mm X 35 cm
LT-206.111	Ø 6 mm X 15 cm
LT-206.112	Ø 6 mm X 20 cm
LT-206.113	Ø 6 mm X 25 cm
LT-206.114	Ø 6 mm X 30 cm
LT-206.115	Ø 6 mm X 35 cm

GENERAL SUCTION CANNULAS



Pyramid Holes Cannula

This commonly used design has three ports located 180 degrees from the thumb print. It is a slightly more aggressive cannula than the Standard during liposuction procedures.

60cc Tommey hub	
LT-201.41	Ø 2.5 mm X 15 cm
LT-201.42	Ø 2.5 mm X 20 cm
LT-201.43	Ø 3.0 mm X 15 cm
LT-201.44	Ø 3.0 mm X 20 cm
LT-201.45	Ø 3.0 mm X 25 cm
LT-201.46	Ø 4.0 mm X 25 cm
LT-201.47	Ø 4.0 mm X 30 cm
LT-201.48	Ø 4.0 mm X 35 cm
LT-201.49	Ø 4.5 mm X 30 cm
LT-201.50	Ø 4.5 mm X 35 cm

Luer lock (Stainless Steel)	
LT-205.51	Ø 1.5 mm X 12 cm
LT-205.52	Ø 2.0 mm X 15 cm
LT-205.53	Ø 2.0 mm X 20 cm
LT-205.54	Ø 2.5 mm X 15 cm
LT-205.55	Ø 2.5 mm X 20 cm
LT-205.56	Ø 3.0 mm X 15 cm
LT-205.57	Ø 3.0 mm X 20 cm
LT-205.58	Ø 3.0 mm X 25 cm

Treaded Hub	
LT-202.41	Ø 3.0 mm X 15 cm
LT-202.42	Ø 3.0 mm X 25 cm
LT-202.43	Ø 4.0 mm X 15 cm
LT-202.44	Ø 4.0 mm X 25 cm

Power Cannula	
LT-206.121	Ø 2 mm X 15 cm
LT-206.122	Ø 2 mm X 20 cm
LT-206.123	Ø 2 mm X 25 cm
LT-206.124	Ø 2 mm X 30 cm
LT-206.125	Ø 2 mm X 35 cm
LT-206.126	Ø 3 mm X 15 cm
LT-206.127	Ø 3 mm X 20 cm
LT-206.128	Ø 3 mm X 25 cm
LT-206.129	Ø 3 mm X 30 cm
LT-206.130	Ø 3 mm X 35 cm
LT-206.131	Ø 4 mm X 15 cm
LT-206.132	Ø 4 mm X 20 cm
LT-206.133	Ø 4 mm X 25 cm
LT-206.134	Ø 4 mm X 30 cm
LT-206.135	Ø 4 mm X 35 cm
LT-206.136	Ø 5 mm X 15 cm
LT-206.137	Ø 5 mm X 20 cm
LT-206.138	Ø 5 mm X 25 cm
LT-206.139	Ø 5 mm X 30 cm
LT-206.140	Ø 5 mm X 35 cm
LT-206.141	Ø 6 mm X 15 cm
LT-206.142	Ø 6 mm X 20 cm
LT-206.143	Ø 6 mm X 25 cm
LT-206.144	Ø 6 mm X 30 cm
LT-206.145	Ø 6 mm X 35 cm

Super Luer Lock	
LT-203.41	Ø 1.5 mm X 12 cm
LT-203.42	Ø 2.0 mm X 15 cm
LT-203.43	Ø 2.5 mm X 15 cm
LT-203.44	Ø 2.5 mm X 20 cm
LT-203.45	Ø 3.0 mm X 15 cm
LT-203.46	Ø 3.0 mm X 20 cm
LT-203.47	Ø 3.0 mm X 25 cm

Luer lock	
LT-204.51	Ø 1.5 mm X 12 cm
LT-204.52	Ø 2.0 mm X 15 cm
LT-204.53	Ø 2.5 mm X 15 cm
LT-204.54	Ø 2.5 mm X 20 cm
LT-204.55	Ø 3.0 mm X 15 cm
LT-204.56	Ø 3.0 mm X 20 cm
LT-204.57	Ø 3.0 mm X 25 cm

GENERAL SUCTION CANNULAS



Mercedes

This circumferential three port configuration cell friendly liposuction cannula is effective for multi-plane soft tissue harvesting.

60cc Tommey hub

LT-201.51	Ø 2.5 mm X 15 cm
LT-201.52	Ø 2.5 mm X 20 cm
LT-201.53	Ø 3.0 mm X 15 cm
LT-201.54	Ø 3.0 mm X 20 cm
LT-201.55	Ø 3.0 mm X 25 cm
LT-201.56	Ø 4.0 mm X 25 cm
LT-201.57	Ø 4.0 mm X 30 cm
LT-201.58	Ø 4.0 mm X 35 cm
LT-201.59	Ø 4.5 mm X 30 cm
LT-201.60	Ø 4.5 mm X 35 cm

Luer lock (Stainless Steel)

LT-205.61	Ø 1.5 mm X 12 cm
LT-205.62	Ø 2.0 mm X 15 cm
LT-205.63	Ø 2.0 mm X 20 cm
LT-205.64	Ø 2.5 mm X 15 cm
LT-205.65	Ø 2.5 mm X 20 cm
LT-205.66	Ø 3.0 mm X 15 cm
LT-205.67	Ø 3.0 mm X 20 cm
LT-205.68	Ø 3.0 mm X 25 cm

Treaded Hub

LT-202.51	Ø 3.0 mm X 15 cm
LT-202.52	Ø 3.0 mm X 25 cm
LT-202.53	Ø 4.0 mm X 15 cm
LT-202.54	Ø 4.0 mm X 25 cm

Super Luer Lock

LT-203.51	Ø 1.5 mm X 12 cm
LT-203.52	Ø 2.0 mm X 15 cm
LT-203.53	Ø 2.5 mm X 15 cm
LT-203.54	Ø 2.5 mm X 20 cm
LT-203.55	Ø 3.0 mm X 15 cm
LT-203.56	Ø 3.0 mm X 20 cm
LT-203.57	Ø 3.0 mm X 25 cm

Luer lock

LT-204.61	Ø 1.5 mm X 12 cm
LT-204.62	Ø 2.0 mm X 15 cm
LT-204.63	Ø 2.5 mm X 15 cm
LT-204.64	Ø 2.5 mm X 20 cm
LT-204.65	Ø 3.0 mm X 15 cm
LT-204.66	Ø 3.0 mm X 20 cm
LT-204.67	Ø 3.0 mm X 25 cm

Power Cannula

LT-206.151	Ø 2 mm X 15 cm
LT-206.152	Ø 2 mm X 20 cm
LT-206.153	Ø 2 mm X 25 cm
LT-206.154	Ø 2 mm X 30 cm
LT-206.155	Ø 2 mm X 35 cm
LT-206.156	Ø 3 mm X 15 cm
LT-206.157	Ø 3 mm X 20 cm
LT-206.158	Ø 3 mm X 25 cm
LT-206.159	Ø 3 mm X 30 cm
LT-206.160	Ø 3 mm X 35 cm
LT-206.161	Ø 4 mm X 15 cm
LT-206.162	Ø 4 mm X 20 cm
LT-206.163	Ø 4 mm X 25 cm
LT-206.164	Ø 4 mm X 30 cm
LT-206.165	Ø 4 mm X 35 cm
LT-206.166	Ø 5 mm X 15 cm
LT-206.167	Ø 5 mm X 20 cm
LT-206.168	Ø 5 mm X 25 cm
LT-206.169	Ø 5 mm X 30 cm
LT-206.170	Ø 5 mm X 35 cm
LT-206.171	Ø 6 mm X 15 cm
LT-206.172	Ø 6 mm X 20 cm
LT-206.173	Ø 6 mm X 25 cm
LT-206.174	Ø 6 mm X 30 cm
LT-206.175	Ø 6 mm X 35 cm

GENERAL SUCTION CANNULAS



3 Spiral Port Cannula (Carraway Style)

Spiral three port design on this Cell Friendly liposuction cannula is effective for multi-plane adipose tissue harvesting.

60cc Tommey hub	
LT-201.61	Ø 2.5 mm X 15 cm
LT-201.62	Ø 2.5 mm X 20 cm
LT-201.63	Ø 3.0 mm X 15 cm
LT-201.64	Ø 3.0 mm X 20 cm
LT-201.65	Ø 3.0 mm X 25 cm
LT-201.66	Ø 4.0 mm X 25 cm
LT-201.67	Ø 4.0 mm X 30 cm
LT-201.68	Ø 4.0 mm X 35 cm
LT-201.69	Ø 4.5 mm X 30 cm
LT-201.70	Ø 4.5 mm X 35 cm

Treaded Hub	
LT-202.61	Ø 3.0 mm X 15 cm
LT-202.62	Ø 3.0 mm X 25 cm
LT-202.63	Ø 4.0 mm X 15 cm
LT-202.64	Ø 4.0 mm X 25 cm

Super Luer Lock	
LT-203.61	Ø 1.5 mm X 12 cm
LT-203.62	Ø 2.0 mm X 15 cm
LT-203.63	Ø 2.5 mm X 15 cm
LT-203.64	Ø 2.5 mm X 20 cm
LT-203.65	Ø 3.0 mm X 15 cm
LT-203.66	Ø 3.0 mm X 20 cm
LT-203.67	Ø 3.0 mm X 25 cm

Luer lock	
LT-204.71	Ø 1.5 mm X 12 cm
LT-204.72	Ø 2.0 mm X 15 cm
LT-204.73	Ø 2.5 mm X 15 cm
LT-204.74	Ø 2.5 mm X 20 cm
LT-204.75	Ø 3.0 mm X 15 cm
LT-204.76	Ø 3.0 mm X 20 cm
LT-204.77	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)	
LT-205.71	Ø 1.5 mm X 12 cm
LT-205.72	Ø 2.0 mm X 15 cm
LT-205.73	Ø 2.0 mm X 20 cm
LT-205.74	Ø 2.5 mm X 15 cm
LT-205.75	Ø 2.5 mm X 20 cm
LT-205.76	Ø 3.0 mm X 15 cm
LT-205.77	Ø 3.0 mm X 20 cm
LT-205.78	Ø 3.0 mm X 25 cm

Power Cannula	
LT-206.181	Ø 2 mm X 15 cm
LT-206.182	Ø 2 mm X 20 cm
LT-206.183	Ø 2 mm X 25 cm
LT-206.184	Ø 2 mm X 30 cm
LT-206.185	Ø 2 mm X 35 cm
LT-206.186	Ø 3 mm X 15 cm
LT-206.187	Ø 3 mm X 20 cm
LT-206.188	Ø 3 mm X 25 cm
LT-206.189	Ø 3 mm X 30 cm
LT-206.190	Ø 3 mm X 35 cm
LT-206.191	Ø 4 mm X 15 cm
LT-206.192	Ø 4 mm X 20 cm
LT-206.193	Ø 4 mm X 25 cm
LT-206.194	Ø 4 mm X 30 cm
LT-206.195	Ø 4 mm X 35 cm
LT-206.196	Ø 5 mm X 15 cm
LT-206.197	Ø 5 mm X 20 cm
LT-206.198	Ø 5 mm X 25 cm
LT-206.199	Ø 5 mm X 30 cm
LT-206.200	Ø 5 mm X 35 cm
LT-206.201	Ø 6 mm X 15 cm
LT-206.202	Ø 6 mm X 20 cm
LT-206.203	Ø 6 mm X 25 cm
LT-206.204	Ø 6 mm X 30 cm
LT-206.205	Ø 6 mm X 35 cm

GENERAL SUCTION CANNULAS



Six Spiral Port (Adjust below 3 Spiral Ports)

Six port spiral configuration on this Cell Friendly liposuction cannula harvests maximum tissue volume in minimal time.

60cc smart hub

LT-201.91	Ø 2.5 mm X 15 cm
LT-201.92	Ø 2.5 mm X 20 cm
LT-201.93	Ø 3.0 mm X 15 cm
LT-201.94	Ø 3.0 mm X 25 cm
LT-201.95	Ø 4.0 mm X 25 cm

Luer lock (Stainless Steel)

LT-205.231	Ø 2.5 mm X 15 cm
LT-205.232	Ø 3.0 mm X 15 cm
LT-205.233	Ø 3.0 mm X 25 cm
LT-205.234	Ø 4.0 mm X 15 cm
LT-205.235	Ø 4.0 mm X 25 cm

Treaded Hub

LT-202.81	Ø 3.0 mm X 15 cm
LT-202.82	Ø 3.0 mm X 25 cm
LT-202.83	Ø 4.0 mm X 15 cm
LT-202.84	Ø 4.0 mm X 25 cm

Power Cannula

LT-206.901	2mm OD X 15 cm
LT-206.902	2mm OD X 20 cm
LT-206.903	2mm OD X 25 cm
LT-206.904	2mm OD X 30 cm
LT-206.905	2mm OD X 35 cm
LT-206.906	3mm OD X 15 cm
LT-206.907	3mm OD X 20 cm
LT-206.908	3mm OD X 25 cm
LT-206.909	3mm OD X 30 cm
LT-206.910	3mm OD X 35 cm
LT-206.911	4mm OD X 15 cm
LT-206.912	4mm OD X 20 cm
LT-206.913	4mm OD X 25 cm
LT-206.914	4mm OD X 30 cm
LT-206.915	4mm OD X 35 cm
LT-206.916	5mm OD X 15 cm
LT-206.917	5mm OD X 20 cm
LT-206.918	5mm OD X 25 cm
LT-206.919	5mm OD X 30 cm
LT-206.920	5mm OD X 35 cm
LT-206.921	6mm OD X 15 cm
LT-206.922	6mm OD X 20 cm
LT-206.923	6mm OD X 25 cm
LT-206.924	6mm OD X 30 cm
LT-206.925	6mm OD X 35 cm

Super Luer Lock

LT-203.81	Ø 2.5 mm X 15 cm
LT-203.82	Ø 3.0 mm X 15 cm
LT-203.83	Ø 3.0 mm X 25 cm
LT-203.84	Ø 4.0 mm X 15 cm
LT-203.85	Ø 4.0 mm X 25 cm

Luer lock

LT-204.91	Ø 2.5 mm X 15 cm
LT-204.92	Ø 3.0 mm X 15 cm
LT-204.93	Ø 3.0 mm X 25 cm
LT-204.94	Ø 4.0 mm X 15 cm
LT-204.95	Ø 4.0 mm X 25 cm

GENERAL SUCTION CANNULAS



Three Port Cannula

This body sculpting instrument is excellent for harvesting small aliquots of tissue during liposuction procedures.

60cc Tommey hub

LT-201.71	Ø 2.5 mm X 15 cm
LT-201.72	Ø 2.5 mm X 20 cm
LT-201.73	Ø 3.0 mm X 15 cm
LT-201.74	Ø 3.0 mm X 20 cm
LT-201.75	Ø 3.0 mm X 25 cm
LT-201.76	Ø 4.0 mm X 25 cm
LT-201.77	Ø 4.0 mm X 30 cm
LT-201.78	Ø 4.0 mm X 35 cm
LT-201.79	Ø 4.5 mm X 30 cm
LT-201.80	Ø 4.5 mm X 35 cm

Treaded Hub

LT-202.71	Ø 3.0 mm X 15 cm
LT-202.72	Ø 3.0 mm X 25 cm
LT-202.73	Ø 4.0 mm X 15 cm
LT-202.74	Ø 4.0 mm X 25 cm

Super Luer Lock

LT-203.71	Ø 1.5 mm X 12 cm
LT-203.72	Ø 2.0 mm X 15 cm
LT-203.73	Ø 2.5 mm X 15 cm
LT-203.74	Ø 2.5 mm X 20 cm
LT-203.75	Ø 3.0 mm X 15 cm
LT-203.76	Ø 3.0 mm X 20 cm
LT-203.77	Ø 3.0 mm X 25 cm

Luer lock

LT-204.81	Ø 1.5 mm X 12 cm
LT-204.82	Ø 2.0 mm X 15 cm
LT-204.83	Ø 2.5 mm X 15 cm
LT-204.84	Ø 2.5 mm X 20 cm
LT-204.85	Ø 3.0 mm X 15 cm
LT-204.86	Ø 3.0 mm X 20 cm
LT-204.87	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

LT-205.81	Ø 1.5 mm X 12 cm
LT-205.82	Ø 2.0 mm X 15 cm
LT-205.83	Ø 2.0 mm X 20 cm
LT-205.84	Ø 2.5 mm X 15 cm
LT-205.85	Ø 2.5 mm X 20 cm
LT-205.86	Ø 3.0 mm X 15 cm
LT-205.87	Ø 3.0 mm X 20 cm
LT-205.88	Ø 3.0 mm X 25 cm

Power Cannula

LT-206.211	Ø 2 mm X 15 cm
LT-206.212	Ø 2 mm X 20 cm
LT-206.213	Ø 2 mm X 25 cm
LT-206.214	Ø 2 mm X 30 cm
LT-206.215	Ø 2 mm X 35 cm
LT-206.216	Ø 3 mm X 15 cm
LT-206.217	Ø 3 mm X 20 cm
LT-206.218	Ø 3 mm X 25 cm
LT-206.219	Ø 3 mm X 30 cm
LT-206.220	Ø 3 mm X 35 cm
LT-206.221	Ø 4 mm X 15 cm
LT-206.222	Ø 4 mm X 20 cm
LT-206.223	Ø 4 mm X 25 cm
LT-206.224	Ø 4 mm X 30 cm
LT-206.225	Ø 4 mm X 35 cm
LT-206.226	Ø 5 mm X 15 cm
LT-206.227	Ø 5 mm X 20 cm
LT-206.228	Ø 5 mm X 25 cm
LT-206.229	Ø 5 mm X 30 cm
LT-206.230	Ø 5 mm X 35 cm
LT-206.231	Ø 6 mm X 15 cm
LT-206.232	Ø 6 mm X 20 cm
LT-206.233	Ø 6 mm X 25 cm
LT-206.234	Ø 6 mm X 30 cm
LT-206.235	Ø 6 mm X 35 cm

GENERAL SUCTION CANNULAS



Basket Shaped Cannula

The Cell Friendly Basket Cannula is a versatile liposculpting tool. It can operate as an equalizer without vacuum to treat adipose tissue and skin irregularities in revision liposuction and fat transfer procedures. The Basket Cannula is also adept under vacuum as a fat harvesting cannula. For finishing, this cannula can be used as an injector for equalizing and adding adipose tissue to sculpt smoother contours.

60cc Tommey hub

LT-201.81	Ø 3.0 mm X 15 cm
LT-201.82	Ø 3.0 mm X 20 cm
LT-201.83	Ø 3.0 mm X 25 cm
LT-201.84	Ø 4.0 mm X 25 cm
LT-201.85	Ø 4.0 mm X 30 cm
LT-201.86	Ø 4.0 mm X 35 cm

Hole on Tip-Style I

Power Cannula

LT-206.251	Ø 2 mm X 15 cm	LT-206.260	Ø 3 mm X 35 cm	LT-206.268	Ø 5 mm X 25 cm
LT-206.252	Ø 2 mm X 20 cm	LT-206.261	Ø 4 mm X 15 cm	LT-206.269	Ø 5 mm X 30 cm
LT-206.253	Ø 2 mm X 25 cm	LT-206.262	Ø 4 mm X 20 cm	LT-206.270	Ø 5 mm X 35 cm
LT-206.254	Ø 2 mm X 30 cm	LT-206.263	Ø 4 mm X 25 cm	LT-206.271	Ø 6 mm X 15 cm
LT-206.255	Ø 2 mm X 35 cm	LT-206.264	Ø 4 mm X 30 cm	LT-206.272	Ø 6 mm X 20 cm
LT-206.256	Ø 3 mm X 15 cm	LT-206.265	Ø 4 mm X 35 cm	LT-206.273	Ø 6 mm X 25 cm
LT-206.257	Ø 3 mm X 20 cm	LT-206.266	Ø 5 mm X 15 cm	LT-206.274	Ø 6 mm X 30 cm
LT-206.258	Ø 3 mm X 25 cm	LT-206.267	Ø 5 mm X 20 cm	LT-206.275	Ø 6 mm X 35 cm
LT-206.259	Ø 3 mm X 30 cm				

Hole on tip-Style II

Power Cannula

LT-206.281	Ø 2 mm X 15 cm	LT-206.290	Ø 3 mm X 35 cm	LT-206.298	Ø 5 mm X 25 cm
LT-206.282	Ø 2 mm X 20 cm	LT-206.291	Ø 4 mm X 15 cm	LT-206.299	Ø 5 mm X 30 cm
LT-206.283	Ø 2 mm X 25 cm	LT-206.292	Ø 4 mm X 20 cm	LT-206.300	Ø 5 mm X 35 cm
LT-206.284	Ø 2 mm X 30 cm	LT-206.293	Ø 4 mm X 25 cm	LT-206.301	Ø 6 mm X 15 cm
LT-206.285	Ø 2 mm X 35 cm	LT-206.294	Ø 4 mm X 30 cm	LT-206.302	Ø 6 mm X 20 cm
LT-206.286	Ø 3 mm X 15 cm	LT-206.295	Ø 4 mm X 35 cm	LT-206.303	Ø 6 mm X 25 cm
LT-206.287	Ø 3 mm X 20 cm	LT-206.296	Ø 5 mm X 15 cm	LT-206.304	Ø 6 mm X 30 cm
LT-206.288	Ø 3 mm X 25 cm	LT-206.297	Ø 5 mm X 20 cm	LT-206.305	Ø 6 mm X 35 cm
LT-206.289	Ø 3 mm X 30 cm				

GENERAL SUCTION CANNULAS



Deep Sharp Cut Dissector

Power Cannula					
LT-206.311	Ø 2 mm X 15 cm	LT-206.320	Ø 3 mm X 35 cm	LT-206.328	Ø 5 mm X 25 cm
LT-206.312	Ø 2 mm X 20 cm	LT-206.321	Ø 4 mm X 15 cm	LT-206.329	Ø 5 mm X 30 cm
LT-206.313	Ø 2 mm X 25 cm	LT-206.322	Ø 4 mm X 20 cm	LT-206.330	Ø 5 mm X 35 cm
LT-206.314	Ø 2 mm X 30 cm	LT-206.323	Ø 4 mm X 25 cm	LT-206.331	Ø 6 mm X 15 cm
LT-206.315	Ø 2 mm X 35 cm	LT-206.324	Ø 4 mm X 30 cm	LT-206.332	Ø 6 mm X 20 cm
LT-206.316	Ø 3 mm X 15 cm	LT-206.325	Ø 4 mm X 35 cm	LT-206.333	Ø 6 mm X 25 cm
LT-206.317	Ø 3 mm X 20 cm	LT-206.326	Ø 5 mm X 15 cm	LT-206.334	Ø 6 mm X 30 cm
LT-206.318	Ø 3 mm X 25 cm	LT-206.327	Ø 5 mm X 20 cm	LT-206.335	Ø 6 mm X 35 cm
LT-206.319	Ø 3 mm X 30 cm				



Bi-Port Cannula

Luer lock (Stainless Steel)				
	10cm	15cm	20cm	25cm
2.0 Ø	LT-205.91	LT-205.101	LT-205.111	LT-205.121
2.5 Ø	LT-205.92	LT-205.102	LT-205.112	LT-205.122
3.0 Ø	LT-205.93	LT-205.103	LT-205.113	LT-205.123
3.5 Ø	LT-205.94	LT-205.104	LT-205.114	LT-205.124
4.0 Ø	LT-205.95	LT-205.105	LT-205.115	LT-205.125

Power Cannula					
LT-206.351	Ø 2 mm X 15 cm	LT-206.360	Ø 3 mm X 35 cm	LT-206.368	Ø 5 mm X 25 cm
LT-206.352	Ø 2 mm X 20 cm	LT-206.361	Ø 4 mm X 15 cm	LT-206.369	Ø 5 mm X 30 cm
LT-206.353	Ø 2 mm X 25 cm	LT-206.362	Ø 4 mm X 20 cm	LT-206.370	Ø 5 mm X 35 cm
LT-206.354	Ø 2 mm X 30 cm	LT-206.363	Ø 4 mm X 25 cm	LT-206.371	Ø 6 mm X 15 cm
LT-206.355	Ø 2 mm X 35 cm	LT-206.364	Ø 4 mm X 30 cm	LT-206.372	Ø 6 mm X 20 cm
LT-206.356	Ø 3 mm X 15 cm	LT-206.365	Ø 4 mm X 35 cm	LT-206.373	Ø 6 mm X 25 cm
LT-206.357	Ø 3 mm X 20 cm	LT-206.366	Ø 5 mm X 15 cm	LT-206.374	Ø 6 mm X 30 cm
LT-206.358	Ø 3 mm X 25 cm	LT-206.367	Ø 5 mm X 20 cm	LT-206.375	Ø 6 mm X 35 cm
LT-206.359	Ø 3 mm X 30 cm				

GENERAL SUCTION CANNULAS



Tapered Tip Single Central Hole

Power Cannula

LT-206.391	Ø 2 mm X 15 cm	LT-206.400	Ø 3 mm X 35 cm	LT-206.408	Ø 5 mm X 25 cm
LT-206.392	Ø 2 mm X 20 cm	LT-206.401	Ø 4 mm X 15 cm	LT-206.409	Ø 5 mm X 30 cm
LT-206.393	Ø 2 mm X 25 cm	LT-206.402	Ø 4 mm X 20 cm	LT-206.410	Ø 5 mm X 35 cm
LT-206.394	Ø 2 mm X 30 cm	LT-206.403	Ø 4 mm X 25 cm	LT-206.411	Ø 6 mm X 15 cm
LT-206.395	Ø 2 mm X 35 cm	LT-206.404	Ø 4 mm X 30 cm	LT-206.412	Ø 6 mm X 20 cm
LT-206.396	Ø 3 mm X 15 cm	LT-206.405	Ø 4 mm X 35 cm	LT-206.413	Ø 6 mm X 25 cm
LT-206.397	Ø 3 mm X 20 cm	LT-206.406	Ø 5 mm X 15 cm	LT-206.414	Ø 6 mm X 30 cm
LT-206.398	Ø 3 mm X 25 cm	LT-206.407	Ø 5 mm X 20 cm	LT-206.415	Ø 6 mm X 35 cm
LT-206.399	Ø 3 mm X 30 cm				



Tapered Tip Two Spiral Holes

Power Cannula

LT-206.431	Ø 2 mm X 15 cm	LT-206.440	Ø 3 mm X 35 cm	LT-206.448	Ø 5 mm X 25 cm
LT-206.432	Ø 2 mm X 20 cm	LT-206.441	Ø 4 mm X 15 cm	LT-206.449	Ø 5 mm X 30 cm
LT-206.433	Ø 2 mm X 25 cm	LT-206.442	Ø 4 mm X 20 cm	LT-206.450	Ø 5 mm X 35 cm
LT-206.434	Ø 2 mm X 30 cm	LT-206.443	Ø 4 mm X 25 cm	LT-206.451	Ø 6 mm X 15 cm
LT-206.435	Ø 2 mm X 35 cm	LT-206.444	Ø 4 mm X 30 cm	LT-206.452	Ø 6 mm X 20 cm
LT-206.436	Ø 3 mm X 15 cm	LT-206.445	Ø 4 mm X 35 cm	LT-206.453	Ø 6 mm X 25 cm
LT-206.437	Ø 3 mm X 20 cm	LT-206.446	Ø 5 mm X 15 cm	LT-206.454	Ø 6 mm X 30 cm
LT-206.438	Ø 3 mm X 25 cm	LT-206.447	Ø 5 mm X 20 cm	LT-206.455	Ø 6 mm X 35 cm
LT-206.439	Ø 3 mm X 30 cm				



Fig. 3F

Power Cannula

LT-206.471	Ø 2 mm X 15 cm	LT-206.480	Ø 3 mm X 35 cm	LT-206.488	Ø 5 mm X 25 cm
LT-206.472	Ø 2 mm X 20 cm	LT-206.481	Ø 4 mm X 15 cm	LT-206.489	Ø 5 mm X 30 cm
LT-206.473	Ø 2 mm X 25 cm	LT-206.482	Ø 4 mm X 20 cm	LT-206.490	Ø 5 mm X 35 cm
LT-206.474	Ø 2 mm X 30 cm	LT-206.483	Ø 4 mm X 25 cm	LT-206.491	Ø 6 mm X 15 cm
LT-206.475	Ø 2 mm X 35 cm	LT-206.484	Ø 4 mm X 30 cm	LT-206.492	Ø 6 mm X 20 cm
LT-206.476	Ø 3 mm X 15 cm	LT-206.485	Ø 4 mm X 35 cm	LT-206.493	Ø 6 mm X 25 cm
LT-206.477	Ø 3 mm X 20 cm	LT-206.486	Ø 5 mm X 15 cm	LT-206.494	Ø 6 mm X 30 cm
LT-206.478	Ø 3 mm X 25 cm	LT-206.487	Ø 5 mm X 20 cm	LT-206.495	Ø 6 mm X 35 cm
LT-206.479	Ø 3 mm X 30 cm				

GENERAL SUCTION CANNULAS



Fig. 2Ech

Power Cannula					
LT-206.501	Ø 2 mm X 15 cm	LT-206.510	Ø 3 mm X 35 cm	LT-206.518	Ø 5 mm X 25 cm
LT-206.502	Ø 2 mm X 20 cm	LT-206.511	Ø 4 mm X 15 cm	LT-206.519	Ø 5 mm X 30 cm
LT-206.503	Ø 2 mm X 25 cm	LT-206.512	Ø 4 mm X 20 cm	LT-206.520	Ø 5 mm X 35 cm
LT-206.504	Ø 2 mm X 30 cm	LT-206.513	Ø 4 mm X 25 cm	LT-206.521	Ø 6 mm X 15 cm
LT-206.505	Ø 2 mm X 35 cm	LT-206.514	Ø 4 mm X 30 cm	LT-206.522	Ø 6 mm X 20 cm
LT-206.506	Ø 3 mm X 15 cm	LT-206.515	Ø 4 mm X 35 cm	LT-206.523	Ø 6 mm X 25 cm
LT-206.507	Ø 3 mm X 20 cm	LT-206.516	Ø 5 mm X 15 cm	LT-206.524	Ø 6 mm X 30 cm
LT-206.508	Ø 3 mm X 25 cm	LT-206.517	Ø 5 mm X 20 cm	LT-206.525	Ø 6 mm X 35 cm
LT-206.509	Ø 3 mm X 30 cm				



Two Square-In Line

Luer lock (Stainless Steel)					
	7cm	10cm	15cm	20cm	25cm
2.0 Ø	LT-205.131	LT-205.141	LT-205.151	LT-205.161	LT-205.171
2.5 Ø	LT-205.132	LT-205.142	LT-205.152	LT-205.162	LT-205.172
3.0 Ø	LT-205.133	LT-205.143	LT-205.153	LT-205.163	LT-205.173
3.5 Ø	LT-205.134	LT-205.144	LT-205.154	LT-205.164	LT-205.174
4.0 Ø	LT-205.135	LT-205.145	LT-205.155	LT-205.165	LT-205.175

Power Cannula					
LT-206.551	Ø 2 mm X 15 cm	LT-206.560	Ø 3 mm X 35 cm	LT-206.568	Ø 5 mm X 25 cm
LT-206.552	Ø 2 mm X 20 cm	LT-206.561	Ø 4 mm X 15 cm	LT-206.569	Ø 5 mm X 30 cm
LT-206.553	Ø 2 mm X 25 cm	LT-206.562	Ø 4 mm X 20 cm	LT-206.570	Ø 5 mm X 35 cm
LT-206.554	Ø 2 mm X 30 cm	LT-206.563	Ø 4 mm X 25 cm	LT-206.571	Ø 6 mm X 15 cm
LT-206.555	Ø 2 mm X 35 cm	LT-206.564	Ø 4 mm X 30 cm	LT-206.572	Ø 6 mm X 20 cm
LT-206.556	Ø 3 mm X 15 cm	LT-206.565	Ø 4 mm X 35 cm	LT-206.573	Ø 6 mm X 25 cm
LT-206.557	Ø 3 mm X 20 cm	LT-206.566	Ø 5 mm X 15 cm	LT-206.574	Ø 6 mm X 30 cm
LT-206.558	Ø 3 mm X 25 cm	LT-206.567	Ø 5 mm X 20 cm	LT-206.575	Ø 6 mm X 35 cm
LT-206.559	Ø 3 mm X 30 cm				

GENERAL SUCTION CANNULAS



One Central Port-Two Parallel Ports

Power Cannula					
LT-206.591	Ø 2 mm X 15 cm	LT-206.600	Ø 3 mm X 35 cm	LT-206.608	Ø 5 mm X 25 cm
LT-206.592	Ø 2 mm X 20 cm	LT-206.601	Ø 4 mm X 15 cm	LT-206.609	Ø 5 mm X 30 cm
LT-206.593	Ø 2 mm X 25 cm	LT-206.602	Ø 4 mm X 20 cm	LT-206.610	Ø 5 mm X 35 cm
LT-206.594	Ø 2 mm X 30 cm	LT-206.603	Ø 4 mm X 25 cm	LT-206.611	Ø 6 mm X 15 cm
LT-206.595	Ø 2 mm X 35 cm	LT-206.604	Ø 4 mm X 30 cm	LT-206.612	Ø 6 mm X 20 cm
LT-206.596	Ø 3 mm X 15 cm	LT-206.605	Ø 4 mm X 35 cm	LT-206.613	Ø 6 mm X 25 cm
LT-206.597	Ø 3 mm X 20 cm	LT-206.606	Ø 5 mm X 15 cm	LT-206.614	Ø 6 mm X 30 cm
LT-206.598	Ø 3 mm X 25 cm	LT-206.607	Ø 5 mm X 20 cm	LT-206.615	Ø 6 mm X 35 cm
LT-206.599	Ø 3 mm X 30 cm				



Two Parallel Port- One Central Port

Power Cannula					
LT-206.651	Ø 2 mm X 15 cm	LT-206.660	Ø 3 mm X 35 cm	LT-206.668	Ø 5 mm X 25 cm
LT-206.652	Ø 2 mm X 20 cm	LT-206.661	Ø 4 mm X 15 cm	LT-206.669	Ø 5 mm X 30 cm
LT-206.653	Ø 2 mm X 25 cm	LT-206.662	Ø 4 mm X 20 cm	LT-206.670	Ø 5 mm X 35 cm
LT-206.654	Ø 2 mm X 30 cm	LT-206.663	Ø 4 mm X 25 cm	LT-206.671	Ø 6 mm X 15 cm
LT-206.655	Ø 2 mm X 35 cm	LT-206.664	Ø 4 mm X 30 cm	LT-206.672	Ø 6 mm X 20 cm
LT-206.656	Ø 3 mm X 15 cm	LT-206.665	Ø 4 mm X 35 cm	LT-206.673	Ø 6 mm X 25 cm
LT-206.657	Ø 3 mm X 20 cm	LT-206.666	Ø 5 mm X 15 cm	LT-206.674	Ø 6 mm X 30 cm
LT-206.658	Ø 3 mm X 25 cm	LT-206.667	Ø 5 mm X 20 cm	LT-206.675	Ø 6 mm X 35 cm
LT-206.659	Ø 3 mm X 30 cm				

GENERAL SUCTION CANNULAS



Fig. 3Me

Luer lock (Stainless Steel)					
	7cm	10cm	15cm	20cm	25cm
2.0 Ø	LT-205.181	LT-205.191	LT-205.201	LT-205.211	LT-205.221
2.5 Ø	LT-205.182	LT-205.192	LT-205.202	LT-205.212	LT-205.222
3.0 Ø	LT-205.183	LT-205.193	LT-205.203	LT-205.213	LT-205.223
3.5 Ø	LT-205.184	LT-205.194	LT-205.204	LT-205.214	LT-205.224
4.0 Ø	LT-205.185	LT-205.195	LT-205.205	LT-205.215	LT-205.225

Power Cannula					
LT-206.681	Ø 2 mm X 15 cm	LT-206.690	Ø 3 mm X 35 cm	LT-206.698	Ø 5 mm X 25 cm
LT-206.682	Ø 2 mm X 20 cm	LT-206.691	Ø 4 mm X 15 cm	LT-206.699	Ø 5 mm X 30 cm
LT-206.683	Ø 2 mm X 25 cm	LT-206.692	Ø 4 mm X 20 cm	LT-206.700	Ø 5 mm X 35 cm
LT-206.684	Ø 2 mm X 30 cm	LT-206.693	Ø 4 mm X 25 cm	LT-206.701	Ø 6 mm X 15 cm
LT-206.685	Ø 2 mm X 35 cm	LT-206.694	Ø 4 mm X 30 cm	LT-206.702	Ø 6 mm X 20 cm
LT-206.686	Ø 3 mm X 15 cm	LT-206.695	Ø 4 mm X 35 cm	LT-206.703	Ø 6 mm X 25 cm
LT-206.687	Ø 3 mm X 20 cm	LT-206.696	Ø 5 mm X 15 cm	LT-206.704	Ø 6 mm X 30 cm
LT-206.688	Ø 3 mm X 25 cm	LT-206.697	Ø 5 mm X 20 cm	LT-206.705	Ø 6 mm X 35 cm
LT-206.689	Ø 3 mm X 30 cm				



Fig. 3Hr

Power Cannula					
LT-206.721	Ø 2 mm X 15 cm	LT-206.730	Ø 3 mm X 35 cm	LT-206.738	Ø 5 mm X 25 cm
LT-206.722	Ø 2 mm X 20 cm	LT-206.731	Ø 4 mm X 15 cm	LT-206.739	Ø 5 mm X 30 cm
LT-206.723	Ø 2 mm X 25 cm	LT-206.732	Ø 4 mm X 20 cm	LT-206.740	Ø 5 mm X 35 cm
LT-206.724	Ø 2 mm X 30 cm	LT-206.733	Ø 4 mm X 25 cm	LT-206.741	Ø 6 mm X 15 cm
LT-206.725	Ø 2 mm X 35 cm	LT-206.734	Ø 4 mm X 30 cm	LT-206.742	Ø 6 mm X 20 cm
LT-206.726	Ø 3 mm X 15 cm	LT-206.735	Ø 4 mm X 35 cm	LT-206.743	Ø 6 mm X 25 cm
LT-206.727	Ø 3 mm X 20 cm	LT-206.736	Ø 5 mm X 15 cm	LT-206.744	Ø 6 mm X 30 cm
LT-206.728	Ø 3 mm X 25 cm	LT-206.737	Ø 5 mm X 20 cm	LT-206.745	Ø 6 mm X 35 cm
LT-206.729	Ø 3 mm X 30 cm				



Fig. Riz

Power Cannula					
LT-206.751	Ø 2 mm X 15 cm	LT-206.760	Ø 3 mm X 35 cm	LT-206.768	Ø 5 mm X 25 cm
LT-206.752	Ø 2 mm X 20 cm	LT-206.761	Ø 4 mm X 15 cm	LT-206.769	Ø 5 mm X 30 cm
LT-206.753	Ø 2 mm X 25 cm	LT-206.762	Ø 4 mm X 20 cm	LT-206.770	Ø 5 mm X 35 cm
LT-206.754	Ø 2 mm X 30 cm	LT-206.763	Ø 4 mm X 25 cm	LT-206.771	Ø 6 mm X 15 cm
LT-206.755	Ø 2 mm X 35 cm	LT-206.764	Ø 4 mm X 30 cm	LT-206.772	Ø 6 mm X 20 cm
LT-206.756	Ø 3 mm X 15 cm	LT-206.765	Ø 4 mm X 35 cm	LT-206.773	Ø 6 mm X 25 cm
LT-206.757	Ø 3 mm X 20 cm	LT-206.766	Ø 5 mm X 15 cm	LT-206.774	Ø 6 mm X 30 cm
LT-206.758	Ø 3 mm X 25 cm	LT-206.767	Ø 5 mm X 20 cm	LT-206.775	Ø 6 mm X 35 cm
LT-206.759	Ø 3 mm X 30 cm				

GENERAL SUCTION CANNULAS



Fig. 5Ai

Power Cannula					
LT-206.721	Ø 2 mm X 15 cm	LT-206.730	Ø 3 mm X 35 cm	LT-206.738	Ø 5 mm X 25 cm
LT-206.722	Ø 2 mm X 20 cm	LT-206.731	Ø 4 mm X 15 cm	LT-206.739	Ø 5 mm X 30 cm
LT-206.723	Ø 2 mm X 25 cm	LT-206.732	Ø 4 mm X 20 cm	LT-206.740	Ø 5 mm X 35 cm
LT-206.724	Ø 2 mm X 30 cm	LT-206.733	Ø 4 mm X 25 cm	LT-206.741	Ø 6 mm X 15 cm
LT-206.725	Ø 2 mm X 35 cm	LT-206.734	Ø 4 mm X 30 cm	LT-206.742	Ø 6 mm X 20 cm
LT-206.726	Ø 3 mm X 15 cm	LT-206.735	Ø 4 mm X 35 cm	LT-206.743	Ø 6 mm X 25 cm
LT-206.727	Ø 3 mm X 20 cm	LT-206.736	Ø 5 mm X 15 cm	LT-206.744	Ø 6 mm X 30 cm
LT-206.728	Ø 3 mm X 25 cm	LT-206.737	Ø 5 mm X 20 cm	LT-206.745	Ø 6 mm X 35 cm
LT-206.729	Ø 3 mm X 30 cm				



Fig. Mr

Power Cannula					
LT-206.751	Ø 2 mm X 15 cm	LT-206.760	Ø 3 mm X 35 cm	LT-206.768	Ø 5 mm X 25 cm
LT-206.752	Ø 2 mm X 20 cm	LT-206.761	Ø 4 mm X 15 cm	LT-206.769	Ø 5 mm X 30 cm
LT-206.753	Ø 2 mm X 25 cm	LT-206.762	Ø 4 mm X 20 cm	LT-206.770	Ø 5 mm X 35 cm
LT-206.754	Ø 2 mm X 30 cm	LT-206.763	Ø 4 mm X 25 cm	LT-206.771	Ø 6 mm X 15 cm
LT-206.755	Ø 2 mm X 35 cm	LT-206.764	Ø 4 mm X 30 cm	LT-206.772	Ø 6 mm X 20 cm
LT-206.756	Ø 3 mm X 15 cm	LT-206.765	Ø 4 mm X 35 cm	LT-206.773	Ø 6 mm X 25 cm
LT-206.757	Ø 3 mm X 20 cm	LT-206.766	Ø 5 mm X 15 cm	LT-206.774	Ø 6 mm X 30 cm
LT-206.758	Ø 3 mm X 25 cm	LT-206.767	Ø 5 mm X 20 cm	LT-206.775	Ø 6 mm X 35 cm
LT-206.759	Ø 3 mm X 30 cm				



Fig. P3G

Power Cannula					
LT-206.791	Ø 2 mm X 15 cm	LT-206.800	Ø 3 mm X 35 cm	LT-206.808	Ø 5 mm X 25 cm
LT-206.792	Ø 2 mm X 20 cm	LT-206.801	Ø 4 mm X 15 cm	LT-206.809	Ø 5 mm X 30 cm
LT-206.793	Ø 2 mm X 25 cm	LT-206.802	Ø 4 mm X 20 cm	LT-206.810	Ø 5 mm X 35 cm
LT-206.794	Ø 2 mm X 30 cm	LT-206.803	Ø 4 mm X 25 cm	LT-206.811	Ø 6 mm X 15 cm
LT-206.795	Ø 2 mm X 35 cm	LT-206.804	Ø 4 mm X 30 cm	LT-206.812	Ø 6 mm X 20 cm
LT-206.796	Ø 3 mm X 15 cm	LT-206.805	Ø 4 mm X 35 cm	LT-206.813	Ø 6 mm X 25 cm
LT-206.797	Ø 3 mm X 20 cm	LT-206.806	Ø 5 mm X 15 cm	LT-206.814	Ø 6 mm X 30 cm
LT-206.798	Ø 3 mm X 25 cm	LT-206.807	Ø 5 mm X 20 cm	LT-206.815	Ø 6 mm X 35 cm
LT-206.799	Ø 3 mm X 30 cm				

GENERAL SUCTION CANNULAS




Fig. Py

Power Cannula					
LT-206.831	Ø 2 mm X 15 cm	LT-206.840	Ø 3 mm X 35 cm	LT-206.848	Ø 5 mm X 25 cm
LT-206.832	Ø 2 mm X 20 cm	LT-206.841	Ø 4 mm X 15 cm	LT-206.849	Ø 5 mm X 30 cm
LT-206.833	Ø 2 mm X 25 cm	LT-206.842	Ø 4 mm X 20 cm	LT-206.850	Ø 5 mm X 35 cm
LT-206.834	Ø 2 mm X 30 cm	LT-206.843	Ø 4 mm X 25 cm	LT-206.851	Ø 6 mm X 15 cm
LT-206.835	Ø 2 mm X 35 cm	LT-206.844	Ø 4 mm X 30 cm	LT-206.852	Ø 6 mm X 20 cm
LT-206.836	Ø 3 mm X 15 cm	LT-206.845	Ø 4 mm X 35 cm	LT-206.853	Ø 6 mm X 25 cm
LT-206.837	Ø 3 mm X 20 cm	LT-206.846	Ø 5 mm X 15 cm	LT-206.854	Ø 6 mm X 30 cm
LT-206.838	Ø 3 mm X 25 cm	LT-206.847	Ø 5 mm X 20 cm	LT-206.855	Ø 6 mm X 35 cm
LT-206.839	Ø 3 mm X 30 cm				



Fig. Le

Power Cannula					
LT-206.861	Ø 2 mm X 15 cm	LT-206.870	Ø 3 mm X 35 cm	LT-206.878	Ø 5 mm X 25 cm
LT-206.862	Ø 2 mm X 20 cm	LT-206.871	Ø 4 mm X 15 cm	LT-206.879	Ø 5 mm X 30 cm
LT-206.863	Ø 2 mm X 25 cm	LT-206.872	Ø 4 mm X 20 cm	LT-206.880	Ø 5 mm X 35 cm
LT-206.864	Ø 2 mm X 30 cm	LT-206.873	Ø 4 mm X 25 cm	LT-206.881	Ø 6 mm X 15 cm
LT-206.865	Ø 2 mm X 35 cm	LT-206.874	Ø 4 mm X 30 cm	LT-206.882	Ø 6 mm X 20 cm
LT-206.866	Ø 3 mm X 15 cm	LT-206.875	Ø 4 mm X 35 cm	LT-206.883	Ø 6 mm X 25 cm
LT-206.867	Ø 3 mm X 20 cm	LT-206.876	Ø 5 mm X 15 cm	LT-206.884	Ø 6 mm X 30 cm
LT-206.868	Ø 3 mm X 25 cm	LT-206.877	Ø 5 mm X 20 cm	LT-206.885	Ø 6 mm X 35 cm
LT-206.869	Ø 3 mm X 30 cm				



Harvester Cannulas

HARVESTER CANNULAS



Multiport Speed Harvester- Both End Sharp

A multiport harvester with elevated 1mm microports which yield small enough aliquots of fat tissue to allow for easy and smooth injection with the Tulip 0.7, 0.9 and 1.2mm injectors.

60cc Tommey hub	
LT-301.01	Ø 2.5 mm X 15 cm
LT-301.02	Ø 2.5 mm X 20 cm
LT-301.03	Ø 3.0 mm X 15 cm
LT-301.04	Ø 3.0 mm X 20 cm
LT-301.05	Ø 3.0 mm X 25 cm
LT-301.06	Ø 4.0 mm X 25 cm
LT-301.07	Ø 4.0 mm X 30 cm
LT-301.08	Ø 4.0 mm X 35 cm
LT-301.09	Ø 4.5 mm X 30 cm
LT-301.10	Ø 4.5 mm X 35 cm

Treaded Hub	
LT-302.01	Ø 3.0 mm X 15 cm
LT-302.02	Ø 3.0 mm X 25 cm
LT-302.03	Ø 4.0 mm X 15 cm
LT-302.04	Ø 4.0 mm X 25 cm

Super Luer Lock	
LT-303.01	Ø 2.0 mm X 15 cm
LT-303.02	Ø 2.5 mm X 15 cm
LT-303.03	Ø 2.5 mm X 20 cm
LT-303.04	Ø 3.0 mm X 15 cm
LT-303.05	Ø 3.0 mm X 20 cm
LT-303.06	Ø 3.0 mm X 25 cm

Luer lock	
LT-304.01	Ø 2.0 mm X 15 cm
LT-304.02	Ø 2.5 mm X 15 cm
LT-304.03	Ø 2.5 mm X 20 cm
LT-304.04	Ø 3.0 mm X 15 cm
LT-304.05	Ø 3.0 mm X 20 cm
LT-304.06	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)	
LT-305.01	Ø 2.0 mm X 15 cm
LT-305.02	Ø 2.5 mm X 15 cm
LT-305.03	Ø 2.5 mm X 20 cm
LT-305.04	Ø 3.0 mm X 15 cm
LT-305.05	Ø 3.0 mm X 20 cm
LT-305.06	Ø 3.0 mm X 25 cm

Power Cannula	
LT-306.01	Ø 2 mm X 15 cm
LT-306.02	Ø 2 mm X 20 cm
LT-306.03	Ø 2 mm X 25 cm
LT-306.04	Ø 2 mm X 30 cm
LT-306.05	Ø 2 mm X 35 cm
LT-306.06	Ø 3 mm X 15 cm
LT-306.07	Ø 3 mm X 20 cm
LT-306.08	Ø 3 mm X 25 cm
LT-306.09	Ø 3 mm X 30 cm
LT-306.10	Ø 3 mm X 35 cm
LT-306.11	Ø 4 mm X 15 cm
LT-306.12	Ø 4 mm X 20 cm
LT-306.13	Ø 4 mm X 25 cm
LT-306.14	Ø 4 mm X 30 cm
LT-306.15	Ø 4 mm X 35 cm
LT-306.16	Ø 5 mm X 15 cm
LT-306.17	Ø 5 mm X 20 cm
LT-306.18	Ø 5 mm X 25 cm
LT-306.19	Ø 5 mm X 30 cm
LT-306.20	Ø 5 mm X 35 cm
LT-306.21	Ø 6 mm X 15 cm
LT-306.22	Ø 6 mm X 20 cm
LT-306.23	Ø 6 mm X 25 cm
LT-306.24	Ø 6 mm X 30 cm
LT-306.25	Ø 6 mm X 35 cm

HARVESTER CANNULAS



Multiport Speed Harvester- One End Sharp

A multiport harvester with elevated 1mm microports which yield small enough aliquots of fat tissue to allow for easy and smooth injection with the Tulip 0.7, 0.9 and 1.2mm injectors.

60cc Tommey hub

LT-301.11	Ø 2.5 mm X 15 cm
LT-301.12	Ø 2.5 mm X 20 cm
LT-301.13	Ø 3.0 mm X 15 cm
LT-301.14	Ø 3.0 mm X 20 cm
LT-301.15	Ø 3.0 mm X 25 cm
LT-301.16	Ø 4.0 mm X 25 cm
LT-301.17	Ø 4.0 mm X 30 cm
LT-301.18	Ø 4.0 mm X 35 cm
LT-301.19	Ø 4.5 mm X 30 cm
LT-301.20	Ø 4.5 mm X 35 cm

Treaded Hub

LT-302.11	Ø 3.0 mm X 15 cm
LT-302.12	Ø 3.0 mm X 25 cm
LT-302.13	Ø 4.0 mm X 15 cm
LT-302.14	Ø 4.0 mm X 25 cm

Super Luer Lock

LT-303.11	Ø 2.0 mm X 15 cm
LT-303.12	Ø 2.5 mm X 15 cm
LT-303.13	Ø 2.5 mm X 20 cm
LT-303.14	Ø 3.0 mm X 15 cm
LT-303.15	Ø 3.0 mm X 20 cm
LT-303.16	Ø 3.0 mm X 25 cm

Luer lock

LT-304.11	Ø 2.0 mm X 15 cm
LT-304.12	Ø 2.5 mm X 15 cm
LT-304.13	Ø 2.5 mm X 20 cm
LT-304.14	Ø 3.0 mm X 15 cm
LT-304.15	Ø 3.0 mm X 20 cm
LT-304.16	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

LT-305.11	Ø 2.0 mm X 15 cm
LT-305.12	Ø 2.5 mm X 15 cm
LT-305.13	Ø 2.5 mm X 20 cm
LT-305.14	Ø 3.0 mm X 15 cm
LT-305.15	Ø 3.0 mm X 20 cm
LT-305.16	Ø 3.0 mm X 25 cm

Power Cannula

LT-306.11	Ø 2 mm X 15 cm
LT-306.12	Ø 2 mm X 20 cm
LT-306.13	Ø 2 mm X 25 cm
LT-306.14	Ø 2 mm X 30 cm
LT-306.15	Ø 2 mm X 35 cm
LT-306.16	Ø 3 mm X 15 cm
LT-306.17	Ø 3 mm X 20 cm
LT-306.18	Ø 3 mm X 25 cm
LT-306.19	Ø 3 mm X 30 cm
LT-306.20	Ø 3 mm X 35 cm
LT-306.21	Ø 4 mm X 15 cm
LT-306.22	Ø 4 mm X 20 cm
LT-306.23	Ø 4 mm X 25 cm
LT-306.24	Ø 4 mm X 30 cm
LT-306.25	Ø 4 mm X 35 cm
LT-306.26	Ø 5 mm X 15 cm
LT-306.27	Ø 5 mm X 20 cm
LT-306.28	Ø 5 mm X 25 cm
LT-306.29	Ø 5 mm X 30 cm
LT-306.30	Ø 5 mm X 35 cm
LT-306.31	Ø 6 mm X 15 cm
LT-306.32	Ø 6 mm X 20 cm
LT-306.33	Ø 6 mm X 25 cm
LT-306.34	Ø 6 mm X 30 cm
LT-306.35	Ø 6 mm X 35 cm



Multi Port Harvester - Both End Bevel

A sophisticated, precision-designed harvester that is truly Patient Friendly. The ports in the Sorensen Harvester are machined with alternating 60 degree cutting edges. This unique design is less traumatic than other multi-port harvesters and delivers a fine-grained, high fat content/purity yield that can easily pass through 1.2mm, 0.9mm and 0.7mm blunt injectors.

60cc Tommey hub	
LT-301.21	Ø 2.5 mm X 15 cm
LT-301.22	Ø 2.5 mm X 20 cm
LT-301.23	Ø 3.0 mm X 15 cm
LT-301.24	Ø 3.0 mm X 20 cm
LT-301.25	Ø 3.0 mm X 25 cm
LT-301.26	Ø 4.0 mm X 25 cm
LT-301.27	Ø 4.0 mm X 30 cm
LT-301.28	Ø 4.0 mm X 35 cm
LT-301.29	Ø 4.5 mm X 30 cm
LT-301.30	Ø 4.5 mm X 35 cm

Treaded Hub	
LT-302.21	Ø 3.0 mm X 15 cm
LT-302.22	Ø 3.0 mm X 25 cm
LT-302.23	Ø 4.0 mm X 15 cm
LT-302.24	Ø 4.0 mm X 25 cm

Super Luer Lock	
LT-303.21	Ø 2.0 mm X 15 cm
LT-303.22	Ø 2.5 mm X 15 cm
LT-303.23	Ø 2.5 mm X 20 cm
LT-303.24	Ø 3.0 mm X 15 cm
LT-303.25	Ø 3.0 mm X 20 cm
LT-303.26	Ø 3.0 mm X 25 cm

Luer lock	
LT-304.21	Ø 2.0 mm X 15 cm
LT-304.22	Ø 2.5 mm X 15 cm
LT-304.23	Ø 2.5 mm X 20 cm
LT-304.24	Ø 3.0 mm X 15 cm
LT-304.25	Ø 3.0 mm X 20 cm
LT-304.26	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)	
LT-305.21	Ø 2.0 mm X 15 cm
LT-305.22	Ø 2.5 mm X 15 cm
LT-305.23	Ø 2.5 mm X 20 cm
LT-305.24	Ø 3.0 mm X 15 cm
LT-305.25	Ø 3.0 mm X 20 cm
LT-305.26	Ø 3.0 mm X 25 cm

Power Cannula	
LT-306.21	Ø 2 mm X 15 cm
LT-306.22	Ø 2 mm X 20 cm
LT-306.23	Ø 2 mm X 25 cm
LT-306.24	Ø 2 mm X 30 cm
LT-306.25	Ø 2 mm X 35 cm
LT-306.26	Ø 3 mm X 15 cm
LT-306.27	Ø 3 mm X 20 cm
LT-306.28	Ø 3 mm X 25 cm
LT-306.29	Ø 3 mm X 30 cm
LT-306.30	Ø 3 mm X 35 cm
LT-306.31	Ø 4 mm X 15 cm
LT-306.32	Ø 4 mm X 20 cm
LT-306.33	Ø 4 mm X 25 cm
LT-306.34	Ø 4 mm X 30 cm
LT-306.35	Ø 4 mm X 35 cm
LT-306.36	Ø 5 mm X 15 cm
LT-306.37	Ø 5 mm X 20 cm
LT-306.38	Ø 5 mm X 25 cm
LT-306.39	Ø 5 mm X 30 cm
LT-306.40	Ø 5 mm X 35 cm
LT-306.41	Ø 6 mm X 15 cm
LT-306.42	Ø 6 mm X 20 cm
LT-306.43	Ø 6 mm X 25 cm
LT-306.44	Ø 6 mm X 30 cm
LT-306.45	Ø 6 mm X 35 cm



Injector Cannulas

(Micro Injectors)

MICRO INJECTORS



Sharp Extractor/Injector



This tip works well for needle extraction of soft tissue and as a soft tissue reinjection needle.

Super Luer Lock

LT-403.01	Ø 1.0 mm X 5 cm
LT-403.02	Ø 1.0 mm X 7 cm
LT-403.03	Ø 2.0 mm X 15 cm
LT-403.04	Ø 2.5 mm X 15 cm
LT-403.05	Ø 3.0 mm X 15 cm

Luer lock

LT-404.01	Ø 1.0 mm X 5 cm
LT-404.02	Ø 1.0 mm X 7 cm
LT-404.03	Ø 2.0 mm X 15 cm
LT-404.04	Ø 2.5 mm X 15 cm
LT-404.05	Ø 3.0 mm X 15 cm

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm
1.1Ø	LT-405.01	LT-405.11	LT-405.21	LT-405.31	LT-405.41
1.4Ø	LT-405.02	LT-405.12	LT-405.22	LT-405.32	LT-405.42
1.8Ø	LT-405.03	LT-405.13	LT-405.23	LT-405.33	LT-405.43
2.0Ø	LT-405.04	LT-405.14	LT-405.24	LT-405.34	LT-405.44

Blunt Extractor/Injector



This tip works well for needle extraction of soft tissue and as a soft tissue reinjection needle.

Super Luer Lock

LT-403.11	Ø 2.0 mm X 15 cm
LT-403.12	Ø 2.5 mm X 15 cm
LT-403.13	Ø 3.0 mm X 15 cm

Luer lock

LT-404.11	Ø 2.0 mm X 15 cm
LT-404.12	Ø 2.5 mm X 15 cm
LT-404.13	Ø 3.0 mm X 15 cm

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm
1.1Ø	LT-405.51	LT-405.61	LT-405.71	LT-405.81	LT-405.91
1.4Ø	LT-405.52	LT-405.62	LT-405.72	LT-405.82	LT-405.92
1.8Ø	LT-405.53	LT-405.63	LT-405.73	LT-405.83	LT-405.93
2.0Ø	LT-405.54	LT-405.64	LT-405.74	LT-405.84	LT-405.94

Blunt Tip Injector



Rounded blunt tip design passes through soft tissue easily without the trauma or potential problems of sharp instruments.

Super Luer Lock

LT-403.21	Ø 1.0 mm X 4 cm
LT-403.22	Ø 1.0 mm X 5 cm
LT-403.23	Ø 1.0 mm X 7 cm
LT-403.24	Curved 12 cm
LT-403.25	Curved 15 cm
LT-403.26	Curved 20 cm

Luer lock

LT-404.21	Ø 1.0 mm X 4 cm
LT-404.22	Ø 1.0 mm X 5 cm
LT-404.23	Ø 1.0 mm X 7 cm
LT-404.24	Curved 12 cm
LT-404.25	Curved 15 cm
LT-404.26	Curved 20 cm

Luer lock (Stainless Steel)

	3cm	5cm	7cm	10cm	12cm
1.1Ø	LT-405.101	LT-405.111	LT-405.121	LT-405.131	LT-405.141
1.4Ø	LT-405.102	LT-405.112	LT-405.122	LT-405.132	LT-405.142
1.8Ø	LT-405.103	LT-405.113	LT-405.123	LT-405.133	LT-405.143
2.0Ø	LT-405.104	LT-405.114	LT-405.124	LT-405.134	LT-405.144

MICRO INJECTORS



 **Fig. 2B**

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	LT-405.201	LT-405.211	LT-405.221	LT-405.231	LT-405.241
1.4Ø	LT-405.202	LT-405.212	LT-405.222	LT-405.232	LT-405.242
1.8Ø	LT-405.203	LT-405.213	LT-405.223	LT-405.233	LT-405.243
2.0Ø	LT-405.204	LT-405.214	LT-405.224	LT-405.234	LT-405.244

 **Fig. 1C**

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	LT-405.251	LT-405.261	LT-405.271	LT-405.281	LT-405.291
1.4Ø	LT-405.252	LT-405.262	LT-405.272	LT-405.282	LT-405.292
1.8Ø	LT-405.253	LT-405.263	LT-405.273	LT-405.283	LT-405.293
2.0Ø	LT-405.254	LT-405.264	LT-405.274	LT-405.284	LT-405.294

 **Fig. 1P**

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	LT-405.301	LT-405.311	LT-405.321	LT-405.331	LT-405.341
1.4Ø	LT-405.302	LT-405.312	LT-405.322	LT-405.332	LT-405.342
1.8Ø	LT-405.303	LT-405.313	LT-405.323	LT-405.333	LT-405.343
2.0Ø	LT-405.304	LT-405.314	LT-405.324	LT-405.334	LT-405.344

 **Fig. 2E**

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	LT-405.351	LT-405.361	LT-405.371	LT-405.381	LT-405.391
1.4Ø	LT-405.352	LT-405.362	LT-405.372	LT-405.382	LT-405.392
1.8Ø	LT-405.353	LT-405.363	LT-405.373	LT-405.383	LT-405.393
2.0Ø	LT-405.354	LT-405.364	LT-405.374	LT-405.384	LT-405.394

MICRO INJECTORS



Coleman Liposuction Cannula Fig. MIC I

Luer lock (Stainless Steel)								
	3cm	5cm	7cm	10cm	12cm	15cm	20cm	3.5cm
1.1Ø	LT-405.401	LT-405.411	LT-405.421	LT-405.431	LT-405.441	LT-405.451	LT-405.461	LT-405.471
1.4Ø	LT-405.402	LT-405.412	LT-405.422	LT-405.432	LT-405.442	LT-405.452	LT-405.462	LT-405.472
1.8Ø	LT-405.403	LT-405.413	LT-405.423	LT-405.433	LT-405.443	LT-405.453	LT-405.463	LT-405.473
2.0Ø	LT-405.404	LT-405.414	LT-405.424	LT-405.434	LT-405.444	LT-405.454	LT-405.464	LT-405.474
2.5Ø	LT-405.405	LT-405.415	LT-405.425	LT-405.435	LT-405.445	LT-405.455	LT-405.465	LT-405.475



Coleman Liposuction Cannula Fig. MIC II

Luer lock (Stainless Steel)							
	3cm	5cm	7cm	10cm	12cm	15cm	20cm
1.1Ø	LT-405.501	LT-405.511	LT-405.521	LT-405.531	LT-405.541	LT-405.551	LT-405.561
1.4Ø	LT-405.502	LT-405.512	LT-405.522	LT-405.532	LT-405.542	LT-405.552	LT-405.562
1.8Ø	LT-405.503	LT-405.513	LT-405.523	LT-405.533	LT-405.543	LT-405.553	LT-405.563
2.0Ø	LT-405.504	LT-405.514	LT-405.524	LT-405.534	LT-405.544	LT-405.554	LT-405.564
2.5Ø	LT-405.505	LT-405.515	LT-405.525	LT-405.535	LT-405.545	LT-405.555	LT-405.565



Coleman Liposuction Cannula Spoon Tip Fig. MIC III

This tip works well for needle extraction of soft tissue and as a soft tissue re injection needle.

Luer lock (Stainless Steel)							
	3cm	5cm	7cm	10cm	12cm	15cm	20cm
1.1Ø	LT-405.601	LT-405.611	LT-405.621	LT-405.631	LT-405.641	LT-405.651	LT-405.661
1.4Ø	LT-405.602	LT-405.612	LT-405.622	LT-405.632	LT-405.642	LT-405.652	LT-405.662
1.8Ø	LT-405.603	LT-405.613	LT-405.623	LT-405.633	LT-405.643	LT-405.653	LT-405.663
2.0Ø	LT-405.604	LT-405.614	LT-405.624	LT-405.634	LT-405.644	LT-405.654	LT-405.664
2.5Ø	LT-405.605	LT-405.615	LT-405.625	LT-405.635	LT-405.645	LT-405.655	LT-405.665



Little Injector

Open and more pointed blunt tip is very effective in dense and fibrous tissue.

Super Luer Lock	
LT-403.51	Ø 1.0 mm X 4 cm
LT-403.52	Ø 1.0 mm X 5 cm
LT-403.53	Ø 1.0 mm X 7 cm

Luer lock	
LT-404.51	Ø 1.0 mm X 4 cm
LT-404.52	Ø 1.0 mm X 5 cm
LT-404.53	Ø 1.0 mm X 7 cm

Luer lock (Stainless Steel)					
	3cm	5cm	7cm	10cm	12cm
1.1Ø	LT-405.701	LT-405.711	LT-405.721	LT-405.731	LT-405.741
1.4Ø	LT-405.702	LT-405.712	LT-405.722	LT-405.732	LT-405.742
1.8Ø	LT-405.703	LT-405.713	LT-405.723	LT-405.733	LT-405.743
2.0Ø	LT-405.704	LT-405.714	LT-405.724	LT-405.734	LT-405.744

MICRO INJECTORS



Fig. 1B

Luer lock (Stainless Steel)							
	3cm	5cm	7cm	10cm	12cm	15cm	20cm
1.1Ø	LT-405.801	LT-405.811	LT-405.821	LT-405.831	LT-405.841	LT-405.851	LT-405.861
1.4Ø	LT-405.802	LT-405.812	LT-405.822	LT-405.832	LT-405.842	LT-405.852	LT-405.862
1.8Ø	LT-405.803	LT-405.813	LT-405.823	LT-405.833	LT-405.843	LT-405.853	LT-405.863
2.0Ø	LT-405.804	LT-405.814	LT-405.824	LT-405.834	LT-405.844	LT-405.854	LT-405.864
2.5Ø	LT-405.805	LT-405.815	LT-405.825	LT-405.835	LT-405.845	LT-405.855	LT-405.865

Curved Injectors



Curved Injectors

New Curved Injectors are designed especially for fat transfer procedures to the body. Their unique shape mimics the anatomical curvature of the female torso and, thus, easily produces the natural-looking results that patients desire. Curved Injection Cannulas ameliorate the fat transfer process by providing a direct and precise route to the injection site and by lessening the time required to correctly sculpt curves.



Speciality Cannulas

Speciality Cannula are useful tools for different to operation are such as face or Neck.

SPECIALITY CANNULAS



V-Dissector

This cannula is specifically designed for use in the superficial compartment of the truncal area. The cutting notch is very useful for undermining in facial procedures and in dense, fibrous tissue.

60cc Tommey hub

LT-501.01	Ø 2.5 mm X 15 cm
LT-501.02	Ø 2.5 mm X 20 cm
LT-501.03	Ø 3.0 mm X 15 cm
LT-501.04	Ø 3.0 mm X 20 cm
LT-501.05	Ø 3.0 mm X 25 cm
LT-501.06	Ø 4.0 mm X 25 cm
LT-501.07	Ø 4.0 mm X 30 cm
LT-501.08	Ø 4.0 mm X 35 cm
LT-501.09	Ø 4.5 mm X 30 cm
LT-501.10	Ø 4.5 mm X 35 cm

Treaded Hub

LT-502.01	Ø 3.0 mm X 15 cm
LT-502.02	Ø 3.0 mm X 25 cm
LT-502.03	Ø 4.0 mm X 15 cm
LT-502.04	Ø 4.0 mm X 25 cm

Super Luer Lock

LT-503.01	Ø 2.0 mm X 15 cm
LT-503.02	Ø 2.5 mm X 15 cm
LT-503.03	Ø 2.5 mm X 20 cm
LT-503.04	Ø 3.0 mm X 15 cm
LT-503.05	Ø 3.0 mm X 20 cm
LT-503.06	Ø 3.0 mm X 25 cm

Luer lock

LT-504.01	Ø 2.0 mm X 15 cm
LT-504.02	Ø 2.5 mm X 15 cm
LT-504.03	Ø 2.5 mm X 20 cm
LT-504.04	Ø 3.0 mm X 15 cm
LT-504.05	Ø 3.0 mm X 20 cm
LT-504.06	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

LT-505.01	Ø 2.0 mm X 15 cm
LT-505.02	Ø 2.5 mm X 15 cm
LT-505.03	Ø 2.5 mm X 20 cm
LT-505.04	Ø 3.0 mm X 15 cm
LT-505.05	Ø 3.0 mm X 20 cm
LT-505.06	Ø 3.0 mm X 25 cm

Power Cannula

LT-506.01	Ø 2 mm X 15 cm
LT-506.02	Ø 2 mm X 20 cm
LT-506.03	Ø 2 mm X 25 cm
LT-506.04	Ø 2 mm X 30 cm
LT-506.05	Ø 2 mm X 35 cm
LT-506.06	Ø 3 mm X 15 cm
LT-506.07	Ø 3 mm X 20 cm
LT-506.08	Ø 3 mm X 25 cm
LT-506.09	Ø 3 mm X 30 cm
LT-506.10	Ø 3 mm X 35 cm
LT-506.11	Ø 4 mm X 15 cm
LT-506.12	Ø 4 mm X 20 cm
LT-506.13	Ø 4 mm X 25 cm
LT-506.14	Ø 4 mm X 30 cm
LT-506.15	Ø 4 mm X 35 cm
LT-506.16	Ø 5 mm X 15 cm
LT-506.17	Ø 5 mm X 20 cm
LT-506.18	Ø 5 mm X 25 cm
LT-506.19	Ø 5 mm X 30 cm
LT-506.20	Ø 5 mm X 35 cm
LT-506.21	Ø 6 mm X 15 cm
LT-506.22	Ø 6 mm X 20 cm
LT-506.23	Ø 6 mm X 25 cm
LT-506.24	Ø 6 mm X 30 cm
LT-506.25	Ø 6 mm X 35 cm

SPECIALITY CANNULAS



Flap

The Flap is primarily used as an undermining dissector, and in smaller diameters, for facial tissue removal. The Flap Cannula in a 2.1 mm diameter is very useful for undermining the nasolabial fold prior to autologous tissue injections, as well as for tissue removal in this anatomic area.

60cc Tommey hub

LT-501.11	Ø 2.5 mm X 15 cm
LT-501.12	Ø 2.5 mm X 20 cm
LT-501.13	Ø 3.0 mm X 15 cm
LT-501.14	Ø 3.0 mm X 20 cm
LT-501.15	Ø 3.0 mm X 25 cm
LT-501.16	Ø 4.0 mm X 25 cm
LT-501.17	Ø 4.0 mm X 30 cm
LT-501.18	Ø 4.0 mm X 35 cm
LT-501.19	Ø 4.5 mm X 30 cm
LT-501.20	Ø 4.5 mm X 35 cm

Treaded Hub

LT-502.11	Ø 3.0 mm X 15 cm
LT-502.12	Ø 3.0 mm X 25 cm
LT-502.13	Ø 4.0 mm X 15 cm
LT-502.14	Ø 4.0 mm X 25 cm

Super Luer Lock

LT-503.11	Ø 2.0 mm X 15 cm
LT-503.12	Ø 2.5 mm X 15 cm
LT-503.13	Ø 2.5 mm X 20 cm
LT-503.14	Ø 3.0 mm X 15 cm
LT-503.15	Ø 3.0 mm X 20 cm
LT-503.16	Ø 3.0 mm X 25 cm

Luer lock

LT-504.11	Ø 2.0 mm X 15 cm
LT-504.12	Ø 2.5 mm X 15 cm
LT-504.13	Ø 2.5 mm X 20 cm
LT-504.14	Ø 3.0 mm X 15 cm
LT-504.15	Ø 3.0 mm X 20 cm
LT-504.16	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

LT-505.11	Ø 2.0 mm X 15 cm
LT-505.12	Ø 2.5 mm X 15 cm
LT-505.13	Ø 2.5 mm X 20 cm
LT-505.14	Ø 3.0 mm X 15 cm
LT-505.15	Ø 3.0 mm X 20 cm
LT-505.16	Ø 3.0 mm X 25 cm

Power Cannula

LT-506.51	Ø 2 mm X 15 cm
LT-506.52	Ø 2 mm X 20 cm
LT-506.53	Ø 2 mm X 25 cm
LT-506.54	Ø 2 mm X 30 cm
LT-506.55	Ø 2 mm X 35 cm
LT-506.56	Ø 3 mm X 15 cm
LT-506.57	Ø 3 mm X 20 cm
LT-506.58	Ø 3 mm X 25 cm
LT-506.59	Ø 3 mm X 30 cm
LT-506.60	Ø 3 mm X 35 cm
LT-506.61	Ø 4 mm X 15 cm
LT-506.62	Ø 4 mm X 20 cm
LT-506.63	Ø 4 mm X 25 cm
LT-506.64	Ø 4 mm X 30 cm
LT-506.65	Ø 4 mm X 35 cm
LT-506.66	Ø 5 mm X 15 cm
LT-506.67	Ø 5 mm X 20 cm
LT-506.68	Ø 5 mm X 25 cm
LT-506.69	Ø 5 mm X 30 cm
LT-506.70	Ø 5 mm X 35 cm
LT-506.71	Ø 6 mm X 15 cm
LT-506.72	Ø 6 mm X 20 cm
LT-506.73	Ø 6 mm X 25 cm
LT-506.74	Ø 6 mm X 30 cm
LT-506.75	Ø 6 mm X 35 cm

SPECIALITY CANNULAS



Field Flap

With its elevated leading edge, this cannula is designed to open and expose soft tissue for more precise suctioning along the mandible and sub-mental areas.

60cc Tommey hub

LT-501.21	Ø 2.5 mm X 15 cm
LT-501.22	Ø 2.5 mm X 20 cm
LT-501.23	Ø 3.0 mm X 15 cm
LT-501.24	Ø 3.0 mm X 20 cm
LT-501.25	Ø 3.0 mm X 25 cm
LT-501.26	Ø 4.0 mm X 25 cm
LT-501.27	Ø 4.0 mm X 30 cm
LT-501.28	Ø 4.0 mm X 35 cm
LT-501.29	Ø 4.5 mm X 30 cm
LT-501.30	Ø 4.5 mm X 35 cm

Treaded Hub

LT-502.21	Ø 3.0 mm X 15 cm
LT-502.22	Ø 3.0 mm X 25 cm
LT-502.23	Ø 4.0 mm X 15 cm
LT-502.24	Ø 4.0 mm X 25 cm

Super Luer Lock

LT-503.21	Ø 2.0 mm X 15 cm
LT-503.22	Ø 2.5 mm X 15 cm
LT-503.23	Ø 2.5 mm X 20 cm
LT-503.24	Ø 3.0 mm X 15 cm
LT-503.25	Ø 3.0 mm X 20 cm
LT-503.26	Ø 3.0 mm X 25 cm

Luer lock

LT-504.21	Ø 2.0 mm X 15 cm
LT-504.22	Ø 2.5 mm X 15 cm
LT-504.23	Ø 2.5 mm X 20 cm
LT-504.24	Ø 3.0 mm X 15 cm
LT-504.25	Ø 3.0 mm X 20 cm
LT-504.26	Ø 3.0 mm X 25 cm

Luer lock (Stainless Steel)

LT-505.21	Ø 2.0 mm X 15 cm
LT-505.22	Ø 2.5 mm X 15 cm
LT-505.23	Ø 2.5 mm X 20 cm
LT-505.24	Ø 3.0 mm X 15 cm
LT-505.25	Ø 3.0 mm X 20 cm
LT-505.26	Ø 3.0 mm X 25 cm

Power Cannula

LT-506.101	Ø 2 mm X 15 cm
LT-506.102	Ø 2 mm X 20 cm
LT-506.103	Ø 2 mm X 25 cm
LT-506.104	Ø 2 mm X 30 cm
LT-506.105	Ø 2 mm X 35 cm
LT-506.106	Ø 3 mm X 15 cm
LT-506.107	Ø 3 mm X 20 cm
LT-506.108	Ø 3 mm X 25 cm
LT-506.109	Ø 3 mm X 30 cm
LT-506.110	Ø 3 mm X 35 cm
LT-506.111	Ø 4 mm X 15 cm
LT-506.112	Ø 4 mm X 20 cm
LT-506.113	Ø 4 mm X 25 cm
LT-506.114	Ø 4 mm X 30 cm
LT-506.115	Ø 4 mm X 35 cm
LT-506.116	Ø 5 mm X 15 cm
LT-506.117	Ø 5 mm X 20 cm
LT-506.118	Ø 5 mm X 25 cm
LT-506.119	Ø 5 mm X 30 cm
LT-506.120	Ø 5 mm X 35 cm
LT-506.121	Ø 6 mm X 15 cm
LT-506.122	Ø 6 mm X 20 cm
LT-506.123	Ø 6 mm X 25 cm
LT-506.124	Ø 6 mm X 30 cm
LT-506.125	Ø 6 mm X 35 cm

SINGLE-USE INSTRUMENTS



Super Luer Lock

The Super Luer Lock is on our Single-use instruments.



Tumescant Infiltrator

This cannula is used specifically for the infusion of tumescant solution and for hydraulic dissection. The configuration of the ports allows even distribution of the fluid.

LT-601.01 Ø 2.0 mm X 12cm
LT-601.02 Ø 2.0 mm X 20cm



Tonnard Harvester

A multiport harvester with 1mm elevated sharp microports. The 20 port configuration allows for harvesting both on the forward and backward strokes.

LT-602.01 Ø 3.0 mm X 20cm
LT-602.02 Ø 2.4 mm X 20cm



Tonnard Harvester

Excellent for harvesting small aliquots of tissues.

LT-603.01 Ø 2.5 mm X 15cm
LT-603.02 Ø 2.0 mm X 15cm



Carraway Harvester

Spiral three port design effective for multi-plane tissue harvesting.

LT-604.01 Ø 2.5 mm X 15cm
LT-604.02 Ø 2.0 mm X 10cm
LT-604.03 Ø 2.0 mm X 15cm
LT-604.04 Ø 1.5 mm X 8cm
LT-604.05 Ø 1.5 mm X 10cm



Injector

Round blunt tip design passes through soft tissue easily without the trauma or potential problems of sharp injectors.

LT-605.01 Ø 1.5 mm X 9cm
LT-605.02 Ø 1.5 mm X 7cm
LT-605.03 Ø 1.0 mm X 5cm
LT-605.04 Ø 1.0 mm X 4cm



Spoon Tip Injector

Blunt injector designed with a more open port for rapid injection of larger volumes - a torso and buttock workhorse.

LT-606.01 Ø 3.0 mm X 20cm
LT-606.02 Ø 2.0 mm X 15cm



Miller Speed Harvester

LT-607.01 Ø 2.0 mm X 15cm

Liposuction Accessories

Handles

Connectors

- a) Luer to Luer
- b) Luer to Handle

Transfer Adapters

- a) Luer to Luer 1.2 mm
Luer to Luer 1.4 mm
Luer to Luer 2.4 mm
- b) Anaerobic Transfer
- c) Micro to Nano

Snap Lock

- a) Standard Snap Locks
60cc, 20cc, 10cc

- b) Snap Lock

Fat Transfer Gun

Decanting Stands

Sterilization Boxed

LIPOSUCTION ACCESSORIES


HANDLES FOR CANNULAS

LT-701.01



Luer Lock Handle

LT-702.01
20cc



Adaptor Handle Style I

LIPOSUCTION ACCESSORIES

HANDLES FOR CANNULAS

LT-704.01

with connection thread for cannulas hose connection



LT-705.01

with connection thread for cannulas hose connection



LT-706.01

with connection thread for cannulas hose connection



LT-707.01

with connection thread for cannulas hose connection



LT-750.01

Trocar

The Trocar is a percutaneous punch that is used to make a eussit dna noitartlfini ot roirp noisicni laitnerefmucric extraction. It is designed to maintain a tight seal while minimizing scarring.



LIPOSUCTION ACCESSORIES

NanoTransfer

LT-751.01

NanoTransfer




LT-752.01

Mesh for NanoTransfer



Anaerobic Transfers

- LT-760.01 Luer 1.2mm
- LT-760.02 Luer 1.4mm
- LT-760.03 Luer 2.4mm



LIPOSUCTION ACCESSORIES

LT-771.01

Luer to Luer



Decant Stands

LT-772.02 60cc to Luer

LT-773.03 Luer to Luer



SnapLoks

These syringe snap locks fit into the lunger flanges of the syringe. When the plunger is pulled out the lock snaps open on the barrel lip thus holding the syringe under vacuum.

LT-901.01	60cc
LT-901.02	50cc
LT-901.03	30cc
LT-901.04	20cc
LT-901.05	10cc



LIPOSUCTION ACCESSORIES



Universal Power Injector

LT-801.01 Universal Fat Injector 10cc & 20cc Two in One

LT-801.02 Fat Injector 60cc



Fat Syringe


LT-950.01
The Power Injector allows the surgeon to re-inject soft tissue through large or small diameter injectors with very little effort and maximum control. The Universal Power Injector will hold a 10cc or a 20cc BD syringe or a BD syringe or a 60cc Monoject syringe.

LIPOSUCTION ACCESSORIES



Sterilization Box

LT-1000.01



Syringe Decanting and Injection

LT-707.01	60 cc
LT-707.02	50 cc
LT-707.03	30 cc
LT-707.04	20 cc
LT-707.05	10 cc



LOTUSMED

INSTRUMENTS



www.lotusmedinstruments.com
Email: lotusmedinstruments@gmail.com
Whatsapp & Call: +92 317 7290854